


LIDCO

Hemodynamic Monitoring

- ✓ Optimise your patients hemodynamic status
- ✓ Non-Invasive, Minimally Invasive, Calibrated, Depth of Anesthesia
- ✓ Parameters include CO, SV, SVV, SVR, PPV



250+

Clinically proven
with 250+ Clinical
Papers

Operating Theatre, ICU, Emergency
Department and Other High Risk Areas

UK Brochure



www.lidco.com

Features



The LiDCOunity monitor is a single platform which combines both the LiDCOplus and LiDCOrapid functions. This provides a single solution to monitoring needs throughout the hospital.

The clinician can choose which mode is most appropriate to the clinical situation. The LiDCOunity can be used non-invasively, minimally invasively with a radial arterial line and can be calibrated with the lithium dilution technique.

It is the only technology available that can be calibrated with a standard radial arterial line and without the need for a central line.

The LiDCOunity uses the PulseCO™ algorithm which converts blood pressure to its constituent parts of flow (CO, SV) and resistance (SVR). The PulseCO™ algorithm is scaled to each patient with either the lithium dilution technique or the nomogram using age, height and weight.

Hemodynamic Monitoring for the entire patient pathway

From the ED to the OR to Critical Care and other High Care departments. LiDCOunity has the flexibility to enable continuity of measurement across patient acuity levels



Emergency Department

LiDCO is used in the emergency department to assist with the early identification of sepsis and the resuscitation of trauma patients

- Evaluate hemodynamic status
- Exclude haemorrhage
- Early hypovolemia diagnosis
- Guide fluid resuscitation
- Early identification of Sepsis
- Guide fluid titration of inotropes

Operating Room

The goal in the OR is to optimise fluid and drug therapy prior and during surgery. Successful hemodynamic monitoring in OR reduces the resources need for high-dependency (ICU)

- Switch seamlessly between Non-Invasive and Minimally Invasive
- Measure Depth of Anesthesia using BIS
- Elected bowel surgery, Aortic aneurism, Vascular surgery

ICU

Recent guidelines published by a Task Force of the ESICM and by the Surviving Sepsis Campaign highlight a need for continuous advanced hemodynamic measurements to guide fluid and drug management

- Can be used with an existing radial arterial line
- Monitoring can be started on admission
- Assess if the patient is fluid responsive
- Start appropriate drug therapy

Other High Risk Areas

LiDCO hemodynamic monitoring systems are used within other high-risk areas. The goal in other high-risk areas is to provide continuous blood pressure monitoring during high-risk procedures such as emergency caesareans

- Maternity
- Cath Lab
- Burns
- Transplant
- Other high dependent areas

LiDCOunity delivers

- ✓ 3 in 1 platform for monitoring in any clinical situation;
- ✓ Non-invasive monitoring, minimally invasive monitoring via the radial arterial line, calibrated values;
- ✓ Very easy to set up and use;
- ✓ Designed to be used to allow for early and prompt monitoring in sepsis;
- ✓ The clinician can see why the blood pressure has changed which helps with important clinical decisions;
- ✓ The LiDCOunity also provides parameters which help to decide when to start and stop giving fluid;
- ✓ One disposable for both arterial line and non-invasive which is cost effective;
- ✓ Depth of anesthesia monitoring with BIS™;
- ✓ Refer to the screen guide tab for further insights into how the flexibility of the displays can help meet your needs.

CONTACT

0207 749 1500

 www.lidco.com



LiDCO[®]unity

Minimally Invasive

- ✓ Plug and play from existing vital signs monitor
- ✓ Arterial line input without needing to change your pressure transducer
- ✓ Validated PulseCO™ algorithm reliably tracks hemodynamic changes in the presence of inotropes and vasoactive drugs
- ✓ Beat-to-beat analysis and display of hemodynamic parameters



LiDCO[®]rapid



Non-Invasive

- ✓ Quick and easy to set-up
- ✓ Real-time continuous non-invasive blood pressure (CNAP™) and hemodynamic parameters
- ✓ Proven to be as effective as an arterial line to monitor fluids when used with the PulseCO™ algorithm
- ✓ Dual finger sensor with automatic finger switching for safer non-invasive use



LiDCO[®]MCNAP

Depth of Anesthesia




- ✓ Integrated into the LiDCO
- ✓ Enables clinicians to titrate anesthesia with its hemodynamic effects
- ✓ Stops over-anesthetising, nor under-anesthetising
- ✓ Stops dramatic falls in blood pressure and flow

LiDCO[®]BIS

Calibrate

- ✓ Continuous real-time measurement with lower risk and high precision
- ✓ Calibrate using LiDCO Lithium technology or another absolute cardiac output measurement value
- ✓ Reduced infection risk with less invasive catheters



LiDCO[®]plus



One Disposable

- ✓ Switch monitoring seamlessly with one disposable Smartcard
- ✓ Smartcard carries key patient information between different LiDCO Monitors to ease set-up and monitoring



Screen Guide



Designed to support your clinical decision making

Long Term Trend

Easy interpretation of trends from the start of monitoring, which can be customised to the parameters you need

Short Term Trend

2 minute window for greater focus on the immediate response to interventions




Numeric data display to assist in recording values for routine clinical charts. The chart display Long term, 2-minute short term trend, event response and preload response screens for LiDCOUnity.

Event Response

Marking and monitoring events like a fluid challenge


Preload Response

Window displays preload response values or volume status indicators of: Pulse Pressure Variation (PPV%) and Stroke Volume Variation (SVV%)



Numeric data display to assist in recording values for routine clinical charts. The chart display displays all absolute and index values.


Touch on any point of the history to review hemodynamic values and review key events



LiDCOplus successful lithium dilution calibration



LiDCOplus key physiology targets and events screen



LiDCOplus hemodynamic target screen, helpful for GDMT

Monitor brain activity with the BIS trend screen

CONTACT

0207 749 1500

www.lidco.com



LiDCO Training Programmes



Royal College of Nursing - Cardiac Output Monitoring

The only Royal College of Nursing accredited study day for Cardiac Output Monitoring.

Who should attend?

- ✓ Critical Care Nurses
- ✓ Professional Development Nurses
- ✓ Junior Doctors

Topics Include:

- ✓ Anatomy and physiology of cardiac output
- ✓ Relevant clinical research
- ✓ Practical, hands-on sessions with simulators
- ✓ Competency based assessment



This one day course has been specifically designed to give our delegates the practical skills and interpretive knowledge underpinned by current theory and research to successfully utilise our monitor in the care setting. Nurses will have the practical skills, confidence, and knowledge to use the equipment along with the tools necessary to facilitate their colleagues' development. Upon successful completion, delegates will receive an RCN certificate and knowledge pack to take back to their units where they can be a great resource for the rest of the team.

Hemodynamic Parameters

Parameter	Equation	Normal Range
Normal Hemodynamic Parameters - Adult		
Arterial Blood Pressure (BP)	Systolic (SBP) Diastolic (DBP)	90 - 140 mmHg 60 - 90 mmHg
Mean Arterial Pressure (MAP)	$SBP + (2 \times DBP)/3$	70 - 105 mmHg
Systolic Pressure Variation (SPV)	$(SP_{max} - SP_{min})^{\nabla}$	<5 mmHg unlikely to be preload responsive >5mmHg likely to be preload responsive
Pulse Pressure Variation (PPV)%	$(PP_{max} - PP_{min}) / [(PP_{max} + PP_{min})/2] \times 100^{\nabla}$	<10% unlikely to be preload responsive >13-15% likely to be preload responsive
Stroke Volume Variation (SVV)%	$(SV_{max} - SV_{min}) / [(SV_{max} + SV_{min})/2] \times 100^{\nabla}$	<10% unlikely to be preload responsive >13-15% likely to be preload responsive
∇ = averaged over 10 sec. of BP data updated every 4 beats		
Right Atrial Pressure (RAP)		2 - 6 mmHg
Right Ventricular Pressure (RVP)	Systolic (RVSP) Diastolic (RVDP)	15 - 25 mmHg 0 - 8 mmHg
Pulmonary Artery Pressure (PAP)	Systolic (PASP) Diastolic (PADP)	15 - 25 mmHg 8 - 15 mmHg
Mean Pulmonary Artery Pressure (MPAP)	$[PASP + (2 \times PADP)]/3$	10 - 20 mmHg
Pulmonary Artery Wedge Pressure (PAWP)		6 - 12 mmHg
Left Atrial Pressure (LAP)		6 - 12 mmHg
Cardiac Output (CO)	$HR \times SV/1000$	4.0 - 8.0 l/min
Cardiac Index (CI)	CO/BSA	2.5 - 4.0 l/min/m ²
Stroke Volume (SV)	$CO/HR \times 1000$	60 - 100 ml/beat
Stroke Volume Index (SVI)	$CI/HR \times 1000$	33 - 47 ml/m ² /beat
Systemic Vascular Resistance (SVR)	$80 \times (MAP - RAP)/CO$	800 - 1200 dynes • sec/cm ⁵
Systemic Vascular Resistance Index (SVRI)	$80 \times MAP - RAP/CI$	1970 - 2390 dynes • sec/cm ⁵ /m ²
Pulmonary Vascular Resistance (PVR)	$80 \times (MPAP - PAWP)/CO$	<250 dynes • sec/cm ⁵
Pulmonary Vascular Resistance Index (PVRI)	$80 \times (MPAP - PAWP)/CI$	255 - 285 dynes • sec/cm ⁵ /m ²
Hemodynamic Parameters - Adult		
Left Ventricular Stroke Work (LVSW)	$SV \times (MAP - PAWP) \times 0.0136$	58 - 104 gm-m/beat
Left Ventricular Stroke Work Index (LVSWI)	$SVI \times (MAP - PAWP) \times 0.0136$	50 - 62 gm-m/m ² /beat
Right Ventricular Stroke Work (RVSW)	$SV \times (MPAP - RAP) \times 0.0136$	8 - 16 gm-m/beat
Right Ventricular Stroke Work Index (RVSWI)	$SVI \times (MPAP - RAP) \times 0.0136$	5 - 10 gm-m/m ² /beat
Coronary Artery Perfusion Pressure (CPP)	Diastolic BP - PAWP	60 - 80 mmHg
Right Ventricular End-Diastolic Volume (RVEDV)	SV/EF	100 - 160 ml
Right Ventricular End-Systolic Volume (RVESV)	EDV - SV	50 - 100 ml
Right Ventricular Ejection Fraction (RVEF)	SV/EDV	40 - 60%
Oxygenation Parameters - Adult		
Partial Pressure of Arterial Oxygen (PaO ₂)		80 - 100 mmHg
Partial Pressure of Arterial CO ₂ (PaCO ₂)		35 - 45 mmHg
Bicarbonate (HCO ₃)		22 - 28 mEq/l
pH		7.38 - 7.42
Arterial Oxygen Saturation (SaO ₂)		95 - 100%
Mixed Venous Saturation (SvO ₂)		60 - 80%
Arterial Oxygen Content (CaO ₂)	$(0.0138 \times Hgb \times SaO_2) + (0.0031 \times PaO_2)$	17 - 20 ml/dl
Venous Oxygen Content (CvO ₂)	$(0.0138 \times Hgb \times SvO_2) + (0.0031 \times PvO_2)$	12 - 15 ml/dl
A-V Oxygen Content Difference (C(a-v)O ₂)	$CaO_2 - CvO_2$	4 - 6 ml/dl
Oxygen Delivery (DO ₂)	$CaO_2 \times CO \times 10$	950 - 1150 ml/min
Oxygen Delivery Index (DO ₂ I)	$CaO_2 \times CI \times 10$	500 - 600 ml/min/m ²
Oxygen Consumption (VO ₂)	$(C(a-v)O_2) \times CO \times 10$	200 - 250 ml/min
Oxygen Consumption Index (VO ₂ I)	$(C(a-v)O_2) \times CI \times 10$	120 - 160 ml/min/m ²
Oxygen Extraction Ratio (O ₂ ER)	$[(CaO_2 - CvO_2)/CaO_2] \times 100$	22 - 30%
Oxygen Extraction Index (O ₂ EI)	$[SaO_2 - SvO_2]/SaO_2 \times 100$	20 - 25%

CONTACT

0207 749 1500

www.lidco.com



Clinical Papers

What Our Customers Say...

References	Reduction In	Average odd or risk ratios (confidence interval)	Number of studies
Ripollésa J, Espinosa A, Martínez-Hurtado M, et al. Intraoperative goal directed hemodynamic therapy in non-cardiac surgery: a systematic review and meta-analysis. <i>Journal of Clinical Anesthesia</i> 2016 Feb; 28: 105-115.	Mortality rate	0.63 (CI: 0.42-0.94)	12
Corcoran T. et al. Perioperative Fluid Management Strategies in Major Surgery: A Stratified Meta-Analysis. <i>Anesthesia -Analgesia</i> 2012; 114(3): 640-651.	Acute kidney injury Pneumonia	0.67 (0.46-0.98) 0.74 (0.57-0.96)	23
Gurgel ST, do Nascimento Jr. P. Maintaining Tissue Perfusion in High-Risk Surgical Patients: A Systematic Review of Randomized Clinical Trials. 2011 International Anesthesia Research Society. DOI: 10.1213/ANE.Ob013e3182055384.	Mortality Organ dysfunction	0.67 (0.55-0.82) 0.62 (0.55-0.70)	32
Aya HD, Cecconi M, Hamilton M, et al. Goal directed therapy in cardiac surgery: a systematic review and meta-analysis. <i>British Journal of Anaesthesia</i> , 2013 Apr;110(4):51D-7.	Postoperative complications Hospital length of stay	0.33 (CI: 0.15-0.73) -2.44 (CI: -4.03 to -0.84)	5
Phan T. Ismail H, Heriot AG, et al. Improving Perioperative Outcomes: Fluid Optimization with the Esophageal Doppler Monitor, a meta-analysis and Review. <i>Journal of the American College of Surgeons</i> , 2008 Dec;207(6):935-41.	Length of stay Postoperative morbidity	-2.34 (CI: -2.91 to -1.77) 0.37 (CI: 0.27-0.50)	9
Arulkumaran N, Corredor C, Hamilton MA, et al. Cardiac complications associated with goal-directed therapy in high-risk surgical patients: a meta-analysis. <i>British Journal of Anaesthesia</i> 2014 Apr;112(4):648-59.	Cardiovascular complications Arrhythmias	0.54 (CI: 0.38-0.76) 0.54 (CI: 0.35-0.85)	22
Cecconi M, Corredor C, Arulkumaran N, et al. Clinical review: Goal-directed therapy-what is the evidence in surgical patients? The effect on different risk groups. <i>Critical Care Medicine</i> 2013, 17:209.	Complications	0.45 (CI: 0.34-Q.60)	32
Dalfino L, Giglio MT, Puntillo F, Marucci M, Brienza N. Haemodynamic goal-directed therapy and postoperative infections: earlier is better. A systematic review and meta-analysis. <i>Critical Care Medicine</i> 2011; 15(3): R154.	Surgical site infection Urinary tract infection Pneumonia	0.58 (0.46-0.74) 0.44 (0.22-0.88) 0.71 (0.55-0.92)	26
Grocott MP, Dushianthan A, Hamilton MA. et al. Perioperative increase in global blood flow to explicit defined goals and outcomes after surgery: a Cochrane systematic review <i>British Journal of Anaesthesia</i> 2013;111(4):535-548.	Acute kidney injury Surgical site infection Respiratory failure Total morbidity rate	0.71 (0.57-0.90) 0.65 (0.50-0.84) 0.51 (0.28-0.93) 0.68 0.58-0.80	31
Srinivasa S, Taylor MH, Sammour T, et al. Oesophageal Doppler-guided fluid administration in colorectal surgery: critical appraisal of published clinical trials. <i>Acta Anaesthesiologica Scandinavica</i> 2011; 55(1): 4-13.	Tissue hypoxia	NA	5
Hamilton MA, Cecconi M, Rhodes A. A systematic review and meta-analysis on the use of preemptive hemodynamic intervention to improve postoperative outcomes in moderate and high risk surgical patients. <i>Anesthesia -Analgesia</i> 2011; 112: 1392-402.	Total morbidity rate	0.44 (0.35-0.55)	29
Brienza N, Giglio MT, Marucci M, et al. Does perioperative hemodynamic optimization protect renal function in surgical patients? A meta-analytic study. <i>Critical Care Medicine</i> 2009;37:2079-90.	Acute kidney injury	0.64 (0.50-0.83)	20
Poeze M, Willem M Greve J, Ramsay G. Meta-analysis of hemodynamic optimization: relationship to methodological quality. <i>Critical Care</i> 2005, 9:R771-R779.	Mortality rate	0.61 (0.46-0.81)	30
Giglio MT, Marucci M, Testini M, et al. Goal-directed haemodynamic therapy and gastrointestinal complications in major surgery: a meta-analysis of randomized controlled trials. <i>British Journal of Anaesthesia</i> ; 2009;103(5):637-646.	Minor gastrointestinal complication Major gastrointestinal complication	0.29 (0.17-0.50) 0.42 (0.27-0.65)	16
Bundgaard-Nielsen M, Holte K, Secher NH, et al. Monitoring of peri-operative fluid administration by individualized goal-directed therapy. <i>Acta Anaesthesiologica Scandinavica</i> 2007 Mar;51(3):331-40.	Hospital length of stay Post-op nausea & vomiting Total morbidity rate	NA	9



We have been using LiDCO products since 2010 within our trust, The device is safe, accurate and easy to use and can be set up within 5 minutes by trained staff. it has the diversity to be used on a conscious, preoperative, perioperative and postoperative patients.

The support received by LiDCO is first class and our registered Nurses have benefitted from RCN accredited study day as well as ad hoc training sessions.

Andy Wright, Medical Equipment Educator
Tameside and Glossop Integrated Care NHS Foundation Trust



The training and support provided by LiDCO has been invaluable for our nursing staff, the trainers have been very knowledgeable and enthusiastic. The nursing team are much more confident to initiate LiDCO and are now happy to act on the results.

Liz Walker, Clinical Education Facilitator
ICCU New Cross Hospital



We use LiDCO plus regularly on septic patients, those with refractory hypotension or those with complex fluid balance management.

We encourage the use of LiDCO to be nurse led on the unit and have invested a great deal of time in training staff on monitoring and calibration. We have received fantastic support from the LiDCO Rep who has been key in delivering training and has been particularly helpful in developing of understanding of troubleshooting. She has been very flexible and understands the fluctuating needs of the unit.

I encourage the use of LiDCO at the bedside as I have found it to be helpful in guiding treatment and informing the clinical picture.

Royal Marsden Hospital



We have been using the LiDCO system on our Critical Care unit for several years now. I find it particularly helpful in managing the haemodynamics of the sickest patients we care for guiding my use of fluids and vasoactive medications. It is especially beneficial in managing patients with septic cardiomyopathy with severe LV impairment who benefit from the addition of inotropes giving me real time, operator independent, cardiac output measurement.

Liz Walker, Clinical Education Facilitator
ICCU New Cross Hospital

CONTACT



0207 749 1500



www.lidco.com

Find Out More
www.lidco.com

LIDCOunity

Contact us to arrange

✓ Product Evaluations ✓ Adoption ✓ Training ✓ Education



UK, Europe & ROW Sales
Tel: +44 (0) 20 7749 1500
Email: info@lidco.com

USA Sales
Toll Free: 877-LIDCO11 (877-543-2611)
Email: usinfo@lidco.com

LIDCO 
www.lidco.com