


PAJUNK®

VascularSono / DeltaLong

Precise placement and radiation-free position control of the CVC



Vascular puncture

MADE IN GERMANY

VascularSono – the perfect complement


Central venous and arterial puncture under ultrasound

The use of ultrasound for the installation of central venous catheters is an evidence-based measure to improve patient well-being and to reduce puncture attempts. With VascularSono, PAJUNK® offers a vascular puncture cannula that has been specially developed for ultrasound placement and is characterised by its excellent echogenic properties. Thanks to the proven Cornerstone Technology, ultrasound waves are reflected very well by both the tip and the cannula shaft, even at a steep insertion angle in the in-plane as well as out-of-plane.




➔ *VascularSono: cannula with excellent echogenicity*

Black marking on the cannula hub makes it easier to check the position of the cannula grinding



Areas of application of the VascularSono cannula out-of-plane and in-plane



Out-of-plane view of the VascularSono in the *vena jugularis interna*



In-plane view of the VascularSono in the *vena subclavia*






Echogenic cannula tip

Facet grinding with two inclination angles




Echogenic cannula shaft

Three-dimensional Cornerstone Reflectors



Optimum sonographic visibility of cannula shaft and tip

- ➔ Reflection on 20 mm length and 360° around the cannula
- ➔ Echogenicity especially at steep insertion angle
- ➔ Clear identification independent of the insertion angle



360° arrangement, evenly offset




Three-dimensional Cornerstone Reflectors

One system – two techniques

DeltaLong – The ECG position control system

The DeltaLong ECG position control system from PAJUNK® is based on a simple basic principle: The intra-atrial ECG potentials are directly derived from the catheter tip. This enables a patient-friendly position control which can be carried out quickly, directly and without exposure to radiation.⁴ PAJUNK® supports two alternative techniques for ECG derivation with DeltaLong.



Method 1:

Injection technique

In this method, the intra-atrial ECG is derived via a physiological saline solution injected through the vein catheter.

Easier handling


The three-way stopcock can be turned 360° and makes handling easier. It can be opened and closed in both directions.

Connection cable
SwitchBox



Improved conductivity

The electrical conductivity of the three-way stopcock has been optimised with the help of gold-plated contact electrodes and the enclosed plug.



Syringe 5 ml, LL

Electrically conductive
three-way stopcock


Prevents formation of bubbles

The stopper, which has a membrane and a Luer Lock Connector, allows the solution to be injected without bubbles and so ensures excellent electrical conductivity.

Method 2:

Seldinger technique

The DeltaLong terminal is connected to the horizontal Seldinger wire.




Connection to SwitchBox
via clamp

The DeltaLong derivation system
is connected to the SwitchBox via the
connection line.

DeltaLong SwitchBox

The flexible connection

The DeltaLong SwitchBox from PAJUNK® offers the right connection and a fast, safe switching option for all common ECG monitors. With just a few simple steps it is possible to create an easy-to-use cable diverter from any standard ECG line – either for diverting a surface electrocardiogram or an intra-atrial ECG via a central venous catheter. The SwitchBox enables “seamless switching” between the two signal sources.



Contact protection
for the connector
of the ECG clamp

Red connecting clamp
for patient electrode

Advantages of the DeltaLong ECG position control system


The ECG potentials derived from the catheter tip are evident. Misalignments can be detected and corrected during placement. The time and cost of radiological position monitoring as well as the radiation exposure of the patient are eliminated.



- ➔ *Faster*
- ➔ *Immediately*
- ➔ *Without radiation exposure*
- ➔ *Accepted by the DGAI (German Society for Anaesthesiology and Intensive Care Medicine)⁴*

Exact position control


An ECG is transferred between the catheter tip and a surface electrode. This changes characteristically when the central venous catheter is pushed into the right atrium (P wave). The significant difference of the P wave between vena cava and atrium is the basis for the localisation of the catheter tip. When the catheter is pushed into the right atrium, the P wave **1** increases and returns to normal when the catheter is retracted **2** so that the correct position can be easily verified.^{5, 6}




VascularSono / DeltaLong

At a glance

VascularSono

Product	Suitable for guidewires up to	Size	Item No.	PU
	0.018 inch	21 G x 35 mm	1187-4F035	25
	0.018 inch	21 G x 70 mm	1187-4F070	25
	0.035 inch	18 G x 40 mm	1187-4K040	25
	0.035 inch	18 G x 70 mm	1187-4K070	25
	0.035 inch	18 G x 100 mm	1187-4K100	25

DeltaLong ECG position control system




DeltaLong kit consisting of:

- ECG connecting cable with UNIVERSAL connector
- Luer-Lock syringe 5 ml; 3-part with silicone plunger
- 3-way stopcock, 360° rotatable, with gold-plated contact electrode and encapsulated plug

Product	Item No.	PU
DeltaLong Set for intra-atrial ECG derivation by injection technique, sterile	1151-01-100	10
DeltaLong Clamp Cable for intra-atrial ECG derivation using Seldinger technique, sterile	1151-01-008	10

DeltaLong SwitchBox



Product	Item No.	PU
DeltaLong SwitchBox	1151-01-000	1

*Studies

- Hansen E., Kutz N., Keyl C., Taeger K. ZVK-Lagekontrolle durch EKG-Ableitung über den Einführdraht, *Anästhesiol Intensivmed. Notfallmed. Schmerzther.* 1wwqw998; 33
- Karaaslan D., Altinisik U., Peker T.T., Nayir E., Ozmen S. External jugular vein catheterization using 'intra-atrial electrocardiogram', *Yonsei Med J.* 2009; 30; 50(2): 222-6
- Kerr R.H., Applegate R.L. Accurate placement of the right atrial air aspiration catheter: a descriptive study and prospective trial of intravascular electrocardiography, *Anesth. Analg.* 2006; 103(2): 435-8
- März P., Postel J., Zierl O. Die Lagekontrolle des Cava-Katheters unter Verwendung der intraatrialen EKG-Ableitung, *Anaesth.* 1984; 33: 123-127
- Weißauer W. Der Cava-Katheter aus medico-legaler Sicht, *Anästhesiol Intensivmed. Notfallmed. Schmerzther.* 1998; 33(2): 117-8
- Wilson R.G., Gaer J.A.R. Right atrial electrocardiography in placement of central venous catheters, *the Lancet* 1988; 462-463
- Schummer W., Sakka S.G., Hüttemann E., Reinhart K., Schummer C., *Ultraschall und Lagekontrolle bei der Anlage zentraler Venenkatheter, Anaesthesist* 2009; 58: 677-685 DOI 10.1007/s00101-009-1569-1
- Edgcombe H., Hocking G. Sonographic identification of needle tip by specialists and novices: a blinded comparison of 5 regional block needles in fresh human cadavers, *Reg. Anesth. Pain Med.* 2010 March-April; 35(2): 207-211

PAJUNK® GmbH
 Medizintechnologie
 Karl-Hall-Strasse 1
 D-78187 Geisingen/Germany
 Phone +49 (0) 77 04/92 91-0
 Telefax +49 (0) 77 04/92 91-6 00
 www.pajunk.com

PAJUNK® Medical Produkte GmbH
 D.A.CH • BeNeLux
 Karl-Hall-Strasse 1
 D-78187 Geisingen/Germany
 Phone +49 (0) 77 04/80 08-0
 Telefax +49 (0) 77 04/80 08-150
 www.pajunk.com