



Ask for the 'purple one'™

NEW

# Anchor\* Tissue Retrieval System™ by CONMED



**AIRSEAL® iFS**  
INSUFFLATION MANAGEMENT  
SYSTEM



**MICROLAP®  
INSTRUMENTS**  
LESS INVASIVE APPROACH



**VCARE®**  
UTERINE MANIPULATOR



To learn more about these and  
other innovative products,  
call 800-448-6506 or  
visit [CONMED.com](http://CONMED.com)

# Anchor\* Tissue Retrieval System™ by CONMED

## SIMPLE

- Intuitive design provides ease of use
- Opens automatically when deployed

## INNOVATIVE

- Unique deployment mechanism offers safety & security
- Proprietary bag material to reduce breaks & leaks

## STRONG

- Rip-stop nylon for exceptional strength & quality
- Provides strong & secure seal for safe tissue retrieval

**TRS-ROBO-8**  
125 ml

**TRS100SB2**  
235 ml


**TRS-ROBO-12**  
300 ml

**TRS-VATS-15**  
1550 ml

**TRS175SB2**  
1550 ml

**TRS190SB2**  
1850 ml

**TRS-TV-25**  
4000 ml





525 French Road, Utica, NY 13502

## Ordering Information

### Anchor\* Tissue Retrieval System™

Introducer/ Trocar Size mm	Mouth Diameter cm / in	Bag Length cm / in	Approx Capacity ml	Qty per Box	Catalog Number
8	5.1 / 2.0	14.0 / 5.5	125	5	TRS-ROB0-8
10	5.6 / 2.2	14.0 / 5.5	235	5	TRS100SB2
12	7.6 / 3.0	14.0 / 5.5	300	5	TRS-ROB0-12
15	14.1 / 5.5	25.4 / 9.9	1550	3	TRS-VATS-15
15	14.1 / 5.5	25.4 / 9.9	1550	3	TRS175SB2
15	15.2 / 6.0	26.5 / 10.4	1850	3	TRS190SB2
25	16.5 / 6.5	38.1 / 15.0	4000	3	TRS-TV-25

Customer Service: 1-800-448-6506    CONMED.com  
International Sales: 1-315-797-8375    info@mail.CONMED.com  
Fax: 1-800-438-3051

\*Anchor is a trademark of the Anchor Products Company, Addison, Illinois.