




**interlock**  
Innovation made in Germany

Cannulas

---

Interlock is a leading supplier of disposable plastic cannulas for topical wound irrigation. Based on our extensive experience, our cannulas are produced in a proprietary injection molding process. Compared to conventional cannulas, our product provides the advantage of superior flexibility: Bending the cannula by more than 90° does not damage the cannula walls. Even multiple twists do not affect the lumen diameter and the cannula functionality remains uncompromised. They are an economical alternative to metal cannulas.

- High value product in comparison to other suppliers
- Large product range: Available in different lengths, diameters, shapes – sterile and non-sterile
- Luer-lock-connectors for multiple adaption possibilities
- 100% quality control: Optical and functional quality control of every single cannula assures a residue free processing



### Interlock disposable cannulas made of plastic - non-sterile

for topical wound irrigation

	REF	Operating length in mm	Total length in mm	Inner Ø in mm	Tip	SU
A	05005	35	45	0,8	Button	200 pcs.
B	05001	60	70	0,8	Button	150 pcs.
C	05003	80	100	1,0	Button	150 pcs.
D	05007	75	90	1,0	Blunt	150 pcs.
E	05009	35	55	0,8	Blunt	200 pcs.

### Interlock disposable button cannulas made of plastic - sterile

for topical wound irrigation

	REF	Operating length in mm	Total length in mm	Inner Ø in mm	PZN	SU
A	05004	35	45	0,8	02435819	200 pcs.
A	05053	35	45	0,8	02497826	50 pcs.
A	05023	35	45	0,8	02497832	21 pcs.
B	05000	60	70	0,8	04954268	150 pcs.
B	05050	60	70	0,8	04954274	50 pcs.
B	05021	60	70	0,8	04954280	21 pcs.
C	05002	80	100	1,0	00424680	150 pcs.
C	05052	80	100	1,0	00425260	50 pcs.
C	05022	80	100	1,0	00426466	21 pcs.

### Interlock disposable cannulas made of plastic - blunt - sterile

for topical wound irrigation

	REF	Operating length in mm	Total length in mm	Inner Ø in mm	PZN	SU
D	05006	75	90	1,0	01673540	150 pcs.
D	05056	75	90	1,0	01674479	50 pcs.
D	05026	75	90	1,0	01674485	21 pcs.
E	05008	35	55	0,8	11161723	200 pcs.
E	05058	35	55	0,8	11161717	50 pcs.
E	05028	35	55	0,8	11161700	21 pcs.

### Interlock disposable button cannulas made of metal - sterile

for topical wound irrigation

	REF	Operating length in mm	Total length in mm	Outer Ø in mm	Inner Ø in mm	SU
F	05013	80	95	0,8	0,50	100 pcs.
G	05011	100	115	1,0	0,60	100 pcs.
H	05047	100	116	1,2	0,65	100 pcs.
I	05038	40	55	1,4	1,00	100 pcs.
J	05046	60	75	1,0	0,65	100 pcs.

### Disposable cannulas made of metal - blunt - sterile

for topical wound irrigation

	REF	Operating length in mm	Total length in mm	Outer Ø in mm	Inner Ø in mm	SU
K	05016	80	95	1,5	1,00	100 pcs.
L	05036	80	95	2,5	2,00	100 pcs.
M	05035	80	95	2,8	2,00	100 pcs.

### Disposable cannulas made of metal - blunt - curved - sterile

for topical wound irrigation

	REF	Operating length in mm	Total length in mm	Outer Ø in mm	Inner Ø in mm	SU
N	05030	100	115	0,6	0,35	100 pcs.
O	05031	100	115	0,8	0,50	100 pcs.
P	05032	80	95	1,0	0,60	100 pcs.
Q	05033	90	105	1,0	0,60	100 pcs.
R	05034	90	105	1,5	1,00	50 pcs.
S	05042	80	95	1,6	1,00	50 pcs.
T	05043	80	95	1,8	1,20	50 pcs.
U	05044	80	95	2,0	1,30	50 pcs.
V	05045	80	95	2,4	1,80	50 pcs.

