




interlock  
Innovation made in Germany

OR &  
Patient care

---

2017





**Interlock Medizintechnik GmbH**

Zum Windpark 1 • 23738 Lensahn • Germany

T: +49 4363 905900 • F: +49 4363 90590590

info@interlockmed.com • www.interlockmed.com



## Contents

	004	•	About Interlock
1	015	•	Cannulae
2	021	•	Bowls and containers
3	033	•	OR and patient care
4	053	•	Hand scrub brushes and brush dispensers
5	057	•	Personal protective equipment
6	081	•	Kirschner wires and accessories
	094	•	Index



**interlock**  
Innovation made in Germany



Interlock Medizintechnik GmbH has its headquarters in the North of Germany. For more than 20 years, Interlock has been developing, producing and distributing products for central sterile services departments (CSSD) at its plant in the municipality of Lensahn, north of Lübeck. From here innovations which always focus on user problems are sent all around the world. The Interlock range has grown continuously over the years and now includes products for endoscopy, surgery and patient care. Today more than 70 employees work for Interlock to ensure that our customers are supplied with the items they require quickly and reliably and that their needs are met.

Interlock Medizintechnik is characterised by...

- Self-regulated production for the uncomplicated manufacture of more than 4,000 high-quality products
- A sophisticated logistics concept to guarantee fast delivery
- Innovative new products, which mean that our range is constantly growing
- Certified quality policy which meets both the applicable standards and our customers' requirements
- International demand for our high-quality products and delivery around the world
- Intensive customer care and order processing by our in-house experts
- Germany-wide customer care and consultation service provided by our nationwide field sales force
- Implementing customer requirements and ideas from the market

## Ordering and product advice

If you require advice or simply wish to place an order, you can contact Interlock from Monday to Friday between 8 a.m. and 5 p.m. as follows:

**Telephone:** +49 4363 905900  
**Fax:** +49 4363 90590590  
**E-mail:** [info@interlockmed.com](mailto:info@interlockmed.com)  
**Website:** [www.interlockmed.com](http://www.interlockmed.com)



### Self-regulated production

- The majority of our products are developed for Interlock in Germany and produced exclusively for Interlock
- In-house production and processing of cannulae and sets under sterile conditions in a state-of-the-art clean room for medical devices
- Own tools for the rapid production of individual products
- In-house assembly
- Long-standing cooperation with reliable partner companies and local suppliers





### Highly efficient logistics

- 48-hour concept ensures rapid product availability for our customers
- Excellent product availability ensured by our 6,000 m<sup>2</sup> warehouse in Lensahn, North Germany.
- State-of-the-art and well-equipped high rack warehouse
- Well-structured warehouse management & separate manual store for rapid turnaround of customer orders
- Reliable cooperation with reputable delivery and shipping companies
- Ideal location right next to the A1 motorway






### **Innovative products**

- Broad and ever-growing range with over 4,000 products for CSSD, surgery and endoscopy
- Exclusively developed products make up the majority of our portfolio
- Special solutions based on customer requirements and suggestions from the market



## Certified quality policy

- Our own certified quality management system in accordance with ISO 13485 as well as conformity to current regulations and standards
- Safety-relevant confirmations and technical data sheets for our products
- Internal quality checks and clearly structured work processes
- Fulfilment of all document requirements within the framework of quality assurance, also for our customers
- Selected partner companies who deliver high-quality products in accordance with international industry standards



## Practical solutions


- Innovative product developments based on customer requirements and suggestions
- Fast implementation of product ideas from the market
- Customised printing and laser engraving for a wide range of products
- Customisable products on request





### Germany-wide presence

- A major name for medical devices and accessories throughout Germany
- Established supplier to many hospitals and medical practices
- Trusted supplier to numerous system providers and dealers
- Presence at numerous trade fairs and congresses in the medical device reprocessing field



### **International presence**

- Worldwide distribution of one of the largest product ranges of CSSD accessories in the world, with more than 3,000 customers in over 40 countries
- Participation at national and international trade fairs and congresses focusing on the healthcare sector
- Catalogue and website in both German and English
- Close relationships with suppliers and manufacturers around the world

### Comprehensive customer service

- Expert advice on the use of our products from experienced medical device consultants
- Highly motivated field sales force for direct support and advice at your premises
- Express delivery of urgent orders on request
- Rapid order handling ensures customer satisfaction





Irrigation cannulae are indispensable medical devices in routine clinical practice and are in particularly high demand in the field of topical wound irrigation. At the same time, it is essential to ensure not only 100% safety but also perfect functioning at all times. Thanks to its many years of experience, Interlock is able to offer you an extensive range of extremely high-quality and yet affordable cannulae made of plastic and stainless steel.

- Disposable plastic cannulae which offer unimpeded flow of liquid even when bent at more than 90°
- Cost-efficient alternatives to reprocessible metal cannulae
- 100% quality control for the customers' safety
- Metal cannulae in different lengths, shapes and diameters
- Cannulae in individual sterile packagings
- Versions with and without olive

Plastic **017**

Metal **019**





A



B



C



D



E

**Interlock disposable cannulae made of plastic – non-sterile packaging**

For topical wound irrigation

	REF	Working length in mm	Total length in mm	Inside dia. in mm	Tip	SU
A	05005	35	45	0.8	Button	200 units
B	05001	60	70	0.8	Button	150 units
C	05003	80	100	1.0	Button	150 units
E	05009	35	55	0.8	Blunt	200 units
D	05007	75	90	1.0	Blunt	150 units

**Interlock disposable button cannulae made of plastic – sterile packaging**

For topical wound irrigation

	REF	Working length in mm	Total length in mm	Inside dia. in mm	German pharmacy central number (PZN)	SU
A	05004	35	45	0.8	02435819	200 units
A	05053	35	45	0.8	02497826	50 units
A	05023	35	45	0.8	02497832	21 units
B	05000	60	70	0.8	04954268	150 units
B	05050	60	70	0.8	04954274	50 units
B	05021	60	70	0.8	04954280	21 units
C	05002	80	100	1.0	00424680	150 units
C	05052	80	100	1.0	00425260	50 units
C	05022	80	100	1.0	00426466	21 units

**Interlock disposable cannulae made of plastic – blunt – sterile packaging**

For topical wound irrigation

	REF	Working length in mm	Total length in mm	Inside dia. in mm	German pharmacy central number (PZN)	SU
E	05008	35	55	0.8	11161723	200 units
E	05058	35	55	0.8	11161717	50 units
E	05028	35	55	0.8	11161700	21 units
D	05006	75	90	1.0	01673540	150 units
D	05056	75	90	1.0	01674479	50 units
D	05026	75	90	1.0	01674485	21 units

A



B



C



D



E



F



G



H



I



J



K



L

M



N



O



P



Q

**Interlock disposable button cannulae made of metal – sterile packaging**

For topical wound irrigation

	REF	Working length in mm	Total length in mm	Outside dia. in mm	Inside dia. in mm	SU
A	05038	40	55	1.40	1.20	100 units
B	05046	60	75	0.90	0.65	100 units
C	05013	80	95	0.80	0.50	100 units
D	05011	100	115	1.00	0.60	100 units
E	05047	100	115	1.20	0.65	100 units

**Interlock disposable cannulae made of metal – blunt – sterile packaging**

For topical wound irrigation

	REF	Working length in mm	Total length in mm	Outside dia. in mm	Inside dia. in mm	SU
F	05016	80	95	1.60	1.00	100 units
G	05036	80	95	2.50	2.00	100 units
H	05035	80	95	2.80	2.00	100 units

**Interlock disposable cannulae made of metal – angled – sterile packaging**

For topical wound irrigation

	REF	Working length in mm	Total length in mm	Outside dia. in mm	Inside dia. in mm	SU
I	05032	80	95	1.00	0.60	100 units
J	05042	80	95	1.60	1.00	50 units
K	05043	80	95	1.80	1.20	50 units
L	05044	80	95	2.00	1.30	50 units
M	05045	80	95	2.40	1.80	50 units
N	05033	90	105	1.00	0.60	100 units
O	05034	90	105	1.50	1.00	50 units
P	05030	100	115	0.60	0.35	100 units
Q	05031	100	115	0.80	0.50	100 units



Bowls are used in hospitals and doctor's surgeries for a range of applications including the preparation of irrigation fluids and disinfectant solutions as well as for the disposal of bandages. Interlock has a wide range of plastic bowls and high-quality stainless steel versions for different purposes on offer in this segment. We even pick and pack complete sets for you in our in-house cleanroom compliant with class 7 as per ISO 14644. Interlock's range also includes a collection of waste containers with a safety lock.

- Bowls made of plastic, stainless steel and aluminium in a range of different sizes
- Graduated bowls for precise filling
- Highest quality standards thanks to optimal material properties
- Sterile packed disposable bowls to allow immediate use
- Complete bowl sets for effective work in the OR
- Waste containers in a range of different sizes and designs



A



B



C



D



E



F



G



H

## Bowls and containers

### Interlock disposable plastic bowls – non-sterile

Can be stacked in each other, with graduation, heat-resistant up to 134°C, for single use only

	REF	Capacity in ml	Colour	SU
A	08240	100	Transparent	100 units
A	08236	100	Blue	100 units
A	08237	100	Yellow	100 units
A	08238	100	Red	100 units
A	08245	100	Red	700 units
A	08239	100	Green	100 units
B	08214	150	Transparent	10 units
B	08247	150	Transparent	700 units
B	08203	150	Blue	10 units
B	08248	150	Blue	700 units
B	08204	150	Yellow	10 units
B	08249	150	Yellow	700 units
B	08205	150	Red	10 units
B	08250	150	Red	700 units
C	08215	250	Transparent	10 units
C	08257	250	Transparent	450 units
C	08206	250	Blue	10 units
C	08252	250	Blue	450 units
C	08207	250	Yellow	10 units
C	08256	250	Yellow	450 units
C	08208	250	Red	10 units
C	08254	250	Red	450 units
D	08209	500	Blue	10 units
D	08255	500	Blue	300 units
E	08210	1000	Blue	10 units
F	08211	3000	Blue	1 unit

### Plastic bowls – non-sterile

Heat-resistant up to 134°C, reusable

	REF	Capacity in ml	Specification	Colour	SU
G	08297	280	-	Blue	80 units
H	08280	5000 (dia. x height: 350 x 130mm)	With graduation, stackable	Blue	1 unit





## Bowls and containers

### Interlock disposable plastic bowls – individual sterile packaging

For single use only

	REF	Capacity in ml	Colour	SU
A	08287	100	Transparent	100 units
A	08288	100	Blue	100 units
A	08289	100	Yellow	100 units
A	08290	100	Red	100 units
A	08291	100	Green	100 units
B	08292	150	Transparent	100 units
B	08293	150	Blue	100 units
B	08294	150	Yellow	100 units
B	08295	150	Red	100 units
C	08218	250	Red	100 units
C	08219	250	Blue	100 units
C	08220	250	Yellow	100 units
D	08227	500	Blue	40 units
E	08228	1000	Blue	35 units
F	08285	3000	Blue	15 units
G	08225	1 x 250 1 x 500	Blue	25 sets
G	08226	2 x 250 2 x 500 1 x 1000	Blue	35 sets
H	08229	1 x 250 1 x spatula 175 mm (white)	Blue	50 sets
H	08234	1 x 1000 1 x spatula 175 mm (white)	Blue	35 sets
H	08235	1 x 500 1 x spatula 175 mm (white)	Blue	40 sets



A



B



C



D



E



F



G

### Kidney dishes made of plastic

	REF	Description	Colour	SU
A	08216	260 x 137 x 45 mm, polypropylene (PP), not resistant to acids and alkalis, heat-resistant up to 134°C, suitable for washer-disinfectors	Blue	5 units
B	08241	700 ml, PET, for single use only	White	100 units
B	08281	700 ml, PET, for single use only, individual sterile packaging	White	35 units

### Aluminium bowls – sterile - for single use only

	REF	Dia. top / bottom in mm	Height in mm	Capacity in ml	Specification	SU
C	08284	62/55	15	32	Individual packaging	200 units
C	08296	80/56	34	110	Pack of 4, sterile packaging	100 x 4 units

### Metal bowls – non-sterile

	REF	Dia. top / bottom in mm	Height in mm	Capacity in ml	Material	Specification	SU
C	08223	62/55	15	32	Aluminium	For single use only	200 units
C	08195	64/51	13	28	Aluminium	For single use only	2100 units
C	08202	80/56	34	110	Aluminium	For single use only	2000 units
C	08196	96/80	25	125	Aluminium	For single use only	2560 units
C	08201	99/70	25	110	Aluminium	For single use only	2000 units
D	08230	-	30	60	Stainless steel	Heat-resistant	2 units
D	08231	-	40	120	Stainless steel	Heat-resistant	2 units
D	08232	-	40	175	Stainless steel	Heat-resistant	2 units
D	08233	-	50	235	Stainless steel	Heat-resistant	2 units

### Interlock spatula made of polypropylene (PP)

Can be sterilised, heat-resistant up to 134°C, for single use only

	REF	Length in mm	Colour	SU
E	08217	160	White	50 units
F	08213	175	Uncoloured	12 units
G	08212	190	Blue	50 units



A



B



C



D



E



F



G



H



I



J

K

L

M

### Griffin and other beakers

Graduated

	REF	Capacity in ml	Designation	Material	SU
A	06004	250	Griffin beaker	Polypropylene	6 units
B	06005	400	Griffin beaker	Polypropylene	6 units
C	06006	500	Griffin beaker	Polypropylene	6 units
D	06018	50	Beaker	Porcelain	1 unit
D	06003	100	Beaker	Porcelain	1 unit
D	06001	250	Beaker	Porcelain	1 unit
D	06002	500	Beaker	Porcelain	1 unit

### Mini plastic cartridges with hinged lid

Heat-resistant up to 134°C, for single use, 40 x 28 x 7 mm

	REF	Colour	SU
E	09889	White	50 units
F	09840	Blue	50 units
G	09841	Green	50 units
H	09843	Yellow	50 units

### Instrument trays

Heat-resistant up to 134°C, suitable for washer-disinfectors

	REF	Format in mm	Material	Colour	SU
I	08488	213 x 113 x 10	Stainless steel	Metal	1 unit
I	08489	235 x 165 x 10	Stainless steel	Metal	1 unit
I	08490	288 x 163 x 10	Stainless steel	Metal	1 unit
J	09832	195 x 135 x 23	Plastic	Red	12 units
K	09833	195 x 135 x 23	Plastic	Yellow	12 units
L	09834	195 x 135 x 23	Plastic	Green	12 units
M	09835	195 x 135 x 23	Plastic	Blue	12 units



A



B



C

**abox® disposal box**

For single use only

	REF	Format in mm	Description	Capacity in l	Colour	SU
A	09879	340 x 403 x 312	Without lid	30	Blue	1 unit
A	09880	340 x 403 x 502	Without lid	50	Black	1 unit
A	09881	-	Lid for REF 09879 and REF 09880	-	Yellow	1 unit

**Disposable sharps containers**

Yellow – with transparent lid and safety seal

	REF	Height in mm	Dia. in mm top	Dia. in mm bottom	Capacity in l	SU
B	26011	150	110	85	1.0	100 units
B	26012	150	140	110	1.5	60 units
B	26013	240	140	120	3.0	30 units

**“PRObox” disposal container**

For single use only, heat-resistant up to 134°C

	REF	Format in mm	Description	SU
C	09826	110 x 90 x 18	For used needles, scalpel blades with safety seal	500 units





**vollständig**

**sauber**

**infektiös**

**unsteril**


**interlock**  
interlock made in Germany

### 3 OR and patient care

Daily clinical routine in the OR and other areas of hospitals requires a wide range of aids and tools – Interlock supplies them. Whether you want to dispose of needles and blades safely or equip your instruments with clamp covers, you’re sure to find the perfect solution on our site. The Interlock range is rounded off by bags for patient’s belongings, sterile light handle covers or sieve tray and repair tags.

- Sieve tray tags in a variety of colours and sizes and with an impressive range of labelling
- Products in sterile packaging for safe use in high-risk areas
- Bags for the storage of patients’ belongings
- Other useful aids for the OR, A&E and wards


Sieve tray tags	035
Repair tags	041
Bags	043
Suture aid booties	045
Miscellaneous	047




A



B



C



D

E

F

Sieve tray tags are used to colour-code sieve trays and containers and to communicate simple information between departments. Interlock is the largest supplier of customised sieve tray tags in the world:

- Customised imprints and colours (e.g., loaner set, non-sterile, cleaned + disinfected, etc.) can be produced quickly and even in small quantities
- Heat- and disinfectant-proof material for use during the entire reprocessing cycle
- Stainless steel racks for sieve tray tags help to keep the workplace well organised

#### Interlock sieve tags made of plastic – blank


For single use only, heat-resistant up to 134°C, suitable for washer-disinfectors

	REF	Format in mm	Colour	SU
A	09634	160 x 35	White	250 units
A	09635	160 x 35	Blue	250 units
A	09636	160 x 35	Green	250 units
A	09637	160 x 35	Red	250 units
A	09638	160 x 35	Yellow	250 units
A	09676	160 x 35	Brown	250 units
A	09677	160 x 35	Black	250 units
B	09244	236 x 50	White	250 units
B	09264	236 x 50	Blue	250 units
B	09265	236 x 50	Green	250 units
B	09266	236 x 50	Red	250 units
B	09267	236 x 50	Yellow	250 units
B	09269	236 x 50	Orange	250 units
C	09566	240 x 50	Blue (one-sided)	250 units
C	09567	240 x 50	Green (one-sided)	250 units
C	09568	240 x 50	Red (one-sided)	250 units
C	09569	240 x 50	Yellow (one-sided)	250 units
D	09268	206 x 110	White	250 units
E	09678	100 x 110	"Fir tree", white	250 units

#### Interlock sieve tray tags made of special paper – blank

For single use only, heat-resistant up to 134°C, not suitable for washer-disinfectors


	REF	Format in mm	Colour	SU
D	09599	206 x 110	White	250 units
F	09751	240 x 50	White	250 units



**Interlock sieve tray tags made of plastic with imprint**

For single use only, heat-resistant up to 134°C, suitable for washer-disinfectors

	REF	Format in mm	Colour	Lettering	SU
	A 09563	240 x 50	Green	"Service"	250 units
	B 09570	240 x 50	Green	"URGENT"	250 units
	C 08309	240 x 50	Red	"MRSA"	250 units
NEW	D 08320	240 x 50	Red	"à réparer"	250 units
	E 09564	240 x 50	Red	"Defective Product"	250 units
	F 09755	240 x 50	Red	"Team-Time-out"	250 units
	G 09752	150 x 35	Red	"Sharpen"	250 units
	H 09584	240 x 50	Yellow	"Urgent"	250 units
NEW	I 08328	240 x 50	Yellow	"incomplet"	250 units
NEW	J 08327	240 x 50	Orange	"en dépôt"	250 units
	K 09753	150 x 35	Orange	"Repair"	250 units
NEW	L 08325	240 x 50	Light blue	"Urgent"	250 units
NEW	M 08319	240 x 50	Violet	"En prêt"	250 units
	N 09565	240 x 50	White	"ORSOS"	250 units
NEW	O 08326	240 x 50	Black	"Réservé"	250 units



**Interlock sieve tray tags made of plastic – on reel**

For single use only, white, core diameter: 76 mm, with tear-off perforation, heat-resistant up to 134°C, suitable for washer-disinfectors

Customised imprint possible as of 100 units

	REF	Format in mm	Description	SU
A	09691	160 x 19	Reel with 5000 units	1 reel
B	09399	254 x 42	Reel with 600 units	1 reel

**Interlock PP sieve tray tags**

Heat-resistant up to 134°C, suitable for washer-disinfectors

	REF	Length in mm	Text field size in mm	Colour	SU
C	09679	200	50 x 30	White	100 units
D	09680	200	50 x 30	Blue	100 units
E	09681	200	50 x 30	Green	100 units
F	09682	200	50 x 30	Red	100 units
G	09683	200	50 x 30	Yellow	100 units

**NEW RFID sieve tray tags made of silicone**

Heat-resistant up to 134°C, integrated RFID chip (ISO-15693)

	REF	Colour	SU
H	08965	Yellow	1 unit
H	08966	Green	1 unit
H	08967	Black	1 unit
H	08968	Blue	1 unit
H	08969	White	1 unit

**Stainless steel rack for sieve tray tags**

	REF	Format in mm	Description	SU	
NEW	I	20441	188 x 100 x 255	3 hooks, 4 self-adhesive feet (12.7 x 1.8 mm)	1 unit
NEW	J	20442	368 x 100 x 255	6 hooks, 4 self-adhesive feet (12.7 x 1.8 mm)	1 unit





A



B



C



D



E



F

**Interlock repair tags**

For single use only

	REF	Colour	Imprint	Description	SU
A	01260	Orange	"reparieren / repair"	Made of plastic	100 units
B	02629	Red	"reparieren / repair"	Made of plastic	100 units
C	01261	Yellow	"à réparer / da ripare"	Made of plastic	100 units
D	01278	Red	"schärfen / sharpen"	Made of plastic	100 units
E	01007	White	"defekt / damaged"	Made of Tyvek, washable at up to 90°C, 160 x 30 mm	1000 units
F	04121	White	-	Made of Tyvek, washable at up to 90°C, 160 x 30 mm	1000 units
F	04132	White	-	Made of Tyvek, washable at up to 90°C, 190 x 30 mm	1000 units
F	04151	White	-	Made of PET, heat-resistant up to 134°C, washable at up to 90°C, 190 x 30 mm	500 units



A



C



D



F



G

**Interlock pouches and bags – made of LDPE**


For single use only

	REF	Format in mm	Type	Description	Imprint	SU
<b>A</b>	03016	585 x 605 x 0.05	Carrying bag	<ul style="list-style-type: none"> <li>• For patient's belongings</li> <li>• With straight flap and resealable adhesive seal</li> <li>• Imprint in 4 languages: German English French Italian</li> </ul>	PATIENTENEIGENTUM / PATIENTS PROPERTY / AFFAIRES PERSONNELLES DU PATIENT / EFFETTI PERSONALI DEL PAZIENTE – Name / name / nom / nome: - Adresse / address / adresse / indirizzo - Zimmer / room / chambre / camera:	250 units
<b>C</b>	03038	200 x 300	Bag with flap	<ul style="list-style-type: none"> <li>• For valuables</li> <li>• With permanent adhesion</li> <li>• With 50 mm flap</li> <li>• Black</li> <li>• Imprint in 4 languages: German English French Italian</li> </ul>	ACHTUNG WERTSACHEN / ATTENTION VALUABLES / ATTENTION OBJETS DE VALEUR / ATTENZIONE OGGETTI DI VALORE	1000 units
<b>D</b>	03021	165 x 240 x 0.045	Bag with flap	<ul style="list-style-type: none"> <li>• For dental prostheses</li> <li>• With permanent adhesion</li> <li>• With 45 mm flap</li> <li>• Yellow</li> <li>• Imprint in 4 languages: German English French Italian</li> </ul>	Achtung: Zahnprothese / Attention: dental prosthesis / Attention: prothèse dentaire / Attenzione: protesi dentaria - NAME / NAME / NOM / COGNOME - ADRESSE / ADDRESS / ADRESSE / INDIRIZZO	500 units
<b>F</b>	09094	220 x 240	Waste bag	<ul style="list-style-type: none"> <li>• With adhesive tape on the rear</li> <li>• White</li> <li>• With 40 mm flap</li> </ul>	-	1000 units


**Interlock suture bag**

For single use only

	REF	Format in mm	Material	Description	SU
<b>G</b>	09093	190 x 105 x 275	Smooth, bleached kraft paper on one side (80 g/m <sup>2</sup> )	<ul style="list-style-type: none"> <li>• With glue which can withstand steam sterilisation</li> <li>• Heat-resistant: up to 140°C</li> <li>• Without imprint</li> </ul>	450 units



A



B



C




D



E



F



G



H



### Suture aid booties

For single use only, sterile packaging

For the protection, marking and identification of sutures – The use of suture aid booties is recommended in a range of surgical fields where thin and delicate suture materials are employed. The suture aid booties protect delicate suture materials from damage. The colour coding of the suture aid booties facilitates the finding of the marked sutures.

	REF	Cap size	Colour	Content	SU
A	42026	Standard	Yellow	5 x 2 pairs	20/box
A	42028	Standard	Yellow	5 x 3 pairs	30/box
A	42020	Standard	Yellow	5 x 5 pairs	50/box
A	42022	Standard	Red	5 x 5 pairs	50/box
A	42023	Standard	Blue	5 x 5 pairs	50/box
A	42024	Standard	White	5 x 5 pairs	50/box
A	42025	Standard	Yellow Red White Blue	5 x 5 pairs	50/box
B	42027	Large	Pink	5 x 2 pairs	20/box
B	42021	Large	Pink	5 x 5 pairs	50/box

### Clamp covers made of silicone – sterile packaging

Length: 125 mm – individual packs of 2

	REF	Size	Designation	Colour	SU
C	42001	1/10"	Mini	Red	40/box
D	42002	2/10"	Midi	White	40/box
E	42003	3/10"	Maxi	Blue	20/box

### Vessel loops – sterile packaging

Length: 45 cm – oval profile – packs of 10 x 2 silicone ribbons

	REF	Format in mm	Colour	Designation	SU
F	41210	1.5 x 1.0	Blue	Mini	20/box
F	41220	1.5 x 1.0	Red	Mini	20/box
F	41230	1.5 x 1.0	White	Mini	20/box
F	41240	1.5 x 1.0	Yellow	Mini	20/box
G	41110	2.5 x 1.0	Blue	Maxi	20/box
G	41120	2.5 x 1.0	Red	Maxi	20/box
G	41130	2.5 x 1.0	White	Maxi	20/box
G	41140	2.5 x 1.0	Yellow	Maxi	20/box
H	41000	5.0 x 1.5	Blue	Super-Maxi	20/box



A

B

C



D



E




F



G



H



I



J

**Interlock MagnoMat – magnetic OR mat**

For the safe resting, transferring and organisation of instruments (made of metal) in the OR. Reduces the risk of sharps injuries.

	REF	Format in mm	Colour	SU
	A 26300	300 x 400	Green	1 unit
NEW	B 26319	300 x 400	Red	1 unit
NEW	C 26318	300 x 400	Yellow	1 unit

**Dustproof bags**

For single use

	REF	Format in mm	Material thickness in µ	Description	SU
	D 03057	90 x 325/35	80	Transparent, not printed, conically tapered, with flap, resealable	100 units

**Safety pins – sterile packaging**

	REF	Length in mm	Blister contents	SU
	E 06007	25	2 units	50 blisters/box
	E 06008	38	2 units	50 blisters/box
	E 06009	50	2 units	50 blisters/box
	F 06019	40	1 unit	250 blisters/box
	F 06020	50	1 unit	250 blisters/box

**Interlock needle counters**

Individually packed, sterile, for single use only

For the safe disposal of used sharp, surgical items such as needles, blades and other small parts.

	REF	Description	SU
	G 26201	Needle counter / sharps container with magnet and transparent lid, 10-fold	100 units
	H 26208	Needle counter / sharps container with magnet and foam inlay, 20-fold	100 units
	I 26209	Needle counter / sharps container with 2 magnets, 20-fold	100 units

**Interlock sharps pads**

Pierce-proof, individual sterile packaging, for single use only

For the safe disposal of used sharp, surgical items such as needles, blades and other small parts.

	REF	Format in mm	Description	SU
	J 26304	145 x 120	The adhesive surfaces hold the items securely in place and seal the pads permanently when folded together.	200 units





A



B



C



D

**Interlock signs**

	REF	Format in mm	Imprint	Material	Packaging	Description	SU
<b>A</b>	04217	270 x 70	Team Time Out	Heat-resistant cardboard	Non-sterile	For single use only	75 units
<b>A</b>	04220	270 x 70	Team Time Out	Heat-resistant cardboard	Sterile	For single use only	50 units
<b>B</b>	20436	140 x 60 x 53	-	Stainless steel (silk-matt glossy brushed)	Non-sterile	Reusable	1 unit

*We will be happy to add your customised lettering with our laser writing technique (1 or 2 lines).  
Please simply enquire when placing your order.*

**Fixing hand and accessories**

Heat-resistant


	REF	Dia. in mm	Material thickness in mm	Designation	Description	SU
<b>C</b>	20600	-	0.98	Fixing hand made of stainless steel	Reusable and bendable, incl. 5 silicone rings	1 unit
<b>D</b>	09611	18	3.5	Silicone ring	Suitable for fixing hand made of stainless steel, heat-resistant	25 units



A




B




C



D



E



F




G



H



I



J

**Light handle covers – sterile packaging**

For single use only

	REF	Version	Quantity / sterile packaging	SU
A	06014	Rigid	1	24/box
B	06015	Flexible	1	100/box
B	06016	Flexible	2	80 x 2/box
B	06017	Flexible	3	80 x 3/box

**Microscope and C-bow covers – individual sterile packaging**

For single use only

	REF	Format in cm	Description	Specification	SU
C	03050	160 x 115	OR microscope cover for Zeiss	For lens ring diameter: 65 mm	1 dispenser with 10 units
D	03051	262 x 115	OR microscope cover for Zeiss	For lens ring diameter: 65 mm	1 dispenser with 10 units
E	03052	300 x 115	OR microscope cover for Zeiss	For lens ring diameter: 65 mm	1 dispenser with 10 units
F	03053	235 x 104	C-bow sterile covers for X-ray devices		1 dispenser with 20 units

**Disposable pillow case made of CPE**

	REF	Format in mm	Colour	SU
G	26009	60 x 60	Blue	250 units

**Identification accessories**

	REF	Format in mm	Designation	Description	SU
H	09141	105 x 58	Bed identification label		100 pads with 100 sheets
I	09210	58 x 110	Bed label	With 2 bevelled corners, metal eyelet with attached rubber cord; blank	750 units

**Boards made of POM synthetic material**

- Protect blades
- With rounded corners
- Suitable for low-temperature sterilisation procedures (not suitable for 134°C steam sterilisation)

	REF	Format in mm	Thickness in mm	Colour	SU
J	09436	24 x 19	3	Uncoloured	1 unit
J	09453	24 x 19	10	Uncoloured	1 unit

interlock



interlock  
Innovation made in Germany

interlock  
Innovation made in Germany

# SURGICAL<sup>®</sup> SCRUB

DRY - WITHOUT DETERGENT

REF 09391



2

STERILE EO

40 Bürsten/Schwämme  
zur Handreinigung im Operationssaal  
40 Trockelschwämme  
für Handreinigung im Operationssaal  
40 Nassschwämme  
für Handreinigung im Operationssaal  
40 Spüllinien  
für Handreinigung im Operationssaal  
40 Spüllinien  
für Handreinigung im Operationssaal

2

STERILE EO



## 4 Hand scrub brushes and brush dispensers

When it comes to patient care, staff members' hands are the primary risk factor in terms of the transmission of nosocomial infections. As such, hand hygiene is one of the most important measures employed in hospitals to avoid these infections. To ensure that not only you but also your patient is always optimally protected, Interlock Medizintechnik offers a select range of hand scrub brushes for routine clinical practice in high-risk areas.

- Standard hand scrubs brushes, which are heat resistant up to 134°C
- Disposable hand scrub brushes for gentle cleaning of hands without skin lesions
- Practical dispenser systems, suitable for all our hand scrub brushes



A



B



C



D



E



F



G



H

## Hand scrub brushes and brush dispensers

### Disposable hand scrub brushes in dispenser

- For single use only
- Individual sterile packaging
- Including Velcro for wall mounting
- With sponge and nail cleaner

	REF	Designation	Disinfectant	SU
A	09391	Interlock Surgical Scrub	Without	40 units
B	09323	Medi scrub	Without	25 units
C	09337	Medi scrub	1.1% iodine solution	25 units
D	09367	Medi scrub	4% chlorhexidine solution	25 units

### Brush dispenser holders

Made of stainless steel

	REF	Description	SU
E	09203	Suitable for Interlock Surgical Scrub dispenser (REF 09391)	1 unit
F	09347	Suitable for Medi Scrub dispenser (REF 09323, 09337, 09367)	1 unit

### Hand scrub brushes

- Non-sterile
- Heat-resistant up to 134°C
- Not suitable for washer-disinfectors
- With strong hollow handle

	REF	Bristle hardness	Description	SU
G	09321	Medium	Not suitable for brush dispenser REF 09329	14 units
H	09357	Soft	Integrated nail brush Not suitable for brush dispenser REF 09329	14 units
H	09327	Soft	Integrated nail brush Suitable for brush dispenser REF 09329	12 units









A



B



C



D



E

Our portfolio includes disposable gowns in a range of versions and materials. For example, lightweight and breathable disposable gowns that are also water-resistant and extremely hard-wearing.

- Coated front and sleeves for the wet area and reprocessing
- Inexpensive disposable visitor gowns for short and less critical use
- Our range includes all standard sizes as well as special XXL gowns

### Interlock gowns

For single use only

	REF	Length in cm	Colour	Material thickness	Description	SU	
	<b>A</b>	09401	136	Yellow	30 g/m <sup>2</sup>	<ul style="list-style-type: none"> <li>- Universal size</li> <li>- Made of non-woven</li> <li>- Non-sterile</li> <li>- Breathable back section</li> <li>- Partially coated (front and sleeves)</li> <li>- Water-resistant</li> <li>- For use in the unclean area of the CSSD</li> <li>- To protect against contaminants</li> <li>- Long sleeves with rubber cuff</li> <li>- Covers entire front and back of the upper body</li> </ul>	100 units
<b>NEW</b>	<b>*</b>	09445	136	Yellow	30 g/m <sup>2</sup>	<ul style="list-style-type: none"> <li>- With cuff</li> <li>- Otherwise as REF 09401</li> </ul>	100 units
	<b>B</b>	09402	136	Orange	15 g/m <sup>2</sup>	<ul style="list-style-type: none"> <li>- Universal size</li> <li>- Made of non-woven</li> <li>- Non-sterile</li> <li>- Breathable back section</li> <li>- Partially coated (front and sleeves)</li> <li>- Water-resistant</li> <li>- For use in the unclean area of the CSSD</li> <li>- To protect against contaminants</li> <li>- Long sleeves with rubber cuff</li> <li>- Covers entire front and back of the upper body</li> </ul>	100 units
	<b>C</b>	09437	130	Blue	72 g/m <sup>2</sup>	<ul style="list-style-type: none"> <li>- Mixed non-woven protective gown</li> <li>- Non-sterile</li> <li>- Breathable</li> <li>- Water-resistant</li> <li>- For use in the unclean area of the CSSD</li> <li>- To protect against contaminants</li> <li>- Long sleeves with cuff</li> <li>- Velcro fastener at neck</li> <li>- Ties at the waist</li> <li>- Covers entire front and back of the upper body</li> </ul>	50 units
	<b>D</b>	09442	150	White	30 g/m <sup>2</sup>	<ul style="list-style-type: none"> <li>- XXL</li> <li>- Non-woven protective gown, non-sterile</li> <li>- For single use</li> <li>- Water-resistant</li> <li>- Breathable</li> <li>- Partially coated (front and sleeves)</li> <li>- Long sleeves with rubber cuff</li> </ul>	100 units
	<b>E</b>	09480	120	Dark green	15 g/m <sup>2</sup>	<ul style="list-style-type: none"> <li>- Visitor gown made of non-woven</li> <li>- Non-sterile</li> <li>- For single use</li> <li>- Not coated</li> <li>- Breathable</li> <li>- Long sleeves with rubber cuff</li> </ul>	100 units



Images to follow soon on our website



A



B



C



D

Interlock's range of protective gloves satisfies the requirements of various working areas. Our nitrile gloves meet high quality standards and have been optimised for the tough demands of the CSSD in terms of ergonomics and performance. Flexibility, resistance and comfort are our key quality criteria, not only for disposable gloves but equally for heat-protective gloves.

- Textured fingertips give exceptional tactile feel
- Nitrile gloves with rolled edge for maximum tear strength and a good fit
- Our disposable gloves come in sizes S to XL and the lengths 24 cm, 30 cm and also 40 cm, for demanding tasks such as the reprocessing of flexible endoscopes
- A large selection of heat-protective gloves round off our portfolio

### Heat protective gloves

	REF	Size	Length in cm	Material	Description	SU
A	09031	9	33	Nitrile coating with non-woven lining	- High heat resistance up to 180°C - Allround protection for the hand and wrist - Extremely breathable - Can be washed at up to 60°C	1 pair
A	09034	10	33	Nitrile coating with non-woven lining	See REF 09031	1 pair
B	09223	9	-	Aramid fibres	- Good mechanical strength at temperatures of up to 100°C - Excellent tactility - High degree of cut protection provided by aramid fibres - Good sweat absorption - High wearing comfort - Can be washed at up to 30°C and in a gentle cycle	1 pair
B	09224	10	-	Aramid fibres	See REF 09223	1 pair
C	09263	-	37	Silicone	- Universal size - Blue - Application range: -60°C to +260°C - Suitable for washer-disinfectors - Latex-free	1 unit
D	09368	-	35	Silicone	- Universal size - Light blue - Application range: -40°C to +240°C - Suitable for washer-disinfectors - Latex-free - Textured surface	1 unit
D	09484	-	43	Silicone	See REF 09368	1 unit



A



B



C

**Scilabub protective gloves made of Nomex®**

- Suitable for hot and cold use
- Application range: -180°C to +260°C
- Textile fabric for good wearing comfort
- Do not shrink or melt
- Hard-wearing
- Good chemical resistance
- Washable
- Red

	REF	Size	Length in cm	SU
A	09883	8.5 - 9	52	1 pair
A	09884	9.5 - 10	52	1 pair
A	09885	10.5 - 11	52	1 pair
B	09886	8.5 - 9	30	1 pair
B	09887	9.5 - 10	30	1 pair
B	09888	10.5 - 11	30	1 pair

**Chemical-protective gloves made of a mixture of neoprene and natural latex**

- Excellent grip and tactile feel
- Not flock-lined
- Material thickness 0.45 mm

	REF	Size	SU
C	09033	6.5 - 7	12 pairs
C	09032	7.5 - 8	12 pairs
C	09030	8.5 - 9	12 pairs





A

B

**Protective gloves made of nitrile**

For single use only

- Powder-free – no risk of contamination of samples/foodstuffs
- Free from latex proteins which can trigger allergies
- In contrast to latex, which can become porous and permeable when overstretched, nitrile is extremely elastic and remains impermeable right up to the tear limit
- Nitrile has an excellent wet grip
- The glove adapts to your hand like a second skin, making it very comfortable to wear. Optimal heat dissipation, no heat build-up. Remains comfortable even when worn for extended periods
- Greater protection than with latex or vinyl: Nitrile is very abrasion- and puncture-resistant
- Resistant to many chemicals and solvents as well as oils and grease

	REF	Size	Length in cm	Material thickness in mm	Colour	SU
A	09362	S	24	0.065	Blue	100 units
A	09360	M	24	0.065	Blue	100 units
A	09361	L	24	0.065	Blue	100 units
B	09364	M	30	0.090	White	100 units
B	09363	L	30	0.090	White	100 units
B	09481	XL	30	0.090	White	100 units



**NEW Nitrile protective gloves – 40 cm long**

- For single use only
- White
- Powder- and latex-free
- High puncture and tear resistance (EN 388)
- High degree of protection against chemicals (EN 374)
- Conforms to category 3, irreversible risks (89/686/EEC)
- Can be worn on both hands
- Roughened finger tips for good grip and safe handling of wet products
- Pleasant tactile feel
- Protective class 3

	REF	Size	Length in cm	Colour	SU
A	09588	S	40	White	100/box
A	09589	M	40	White	100/box
A	09590	L	40	White	100/box
A	09591	XL	40	White	100/box
A	09592	XXL	40	White	100/box

**Nitrile protective gloves – 40 cm long**

- For single use only
- Blue
- Hard-wearing, roughened
- Powder- and latex-free
- Protective class 1

	REF	Size	Length in cm	Colour	SU
B	09472	S	40	Blue	50/box
B	09473	M	40	Blue	50/box
B	09474	L	40	Blue	50/box
B	09475	XL	40	Blue	50/box



A



B

**Cleaning gloves made of PE – 90 cm long**

- For single use only
- Made of polyethylene (soft film)
- Material thickness 30 µ
- Elasticated

	REF	Colour	SU
A	09836	Green	50 units
A	09860	Orange	50 units

**Protective sleeves made of PE – 40 cm long**

- For single use only
- Made of polyethylene (PE)
- Material thickness 20 µ
- Moisture-proof
- Hard-wearing
- Elasticated

	REF	Colour	SU
B	09438	White	100 units



A



B



C



D



E



F



G



H



I



J

## Personal protective equipment • Protective face shields

Safety goggles and visors have an important role to play, particularly when reprocessing and working in the wet area. It is not only contaminated instruments which present a risk but also cleaning chemicals, which can cause severe harm. Accordingly, we attach great importance to high-quality products, such as our full-face visors:

- Lightweight safety goggles and visors ensure maximum comfort when worn for extended periods
- Face visors feature anti-fog properties which guarantee a clear view at all times

### Protective face shields

For single use only

	REF	Designation	Description	Colour	SU
A	09081	Interlock goggles frame	For 2-piece safety goggles	Yellow	10 units
B	09431	Interlock goggles frame	For 2-piece safety goggles	Pink	10 units
C	09082	Interlock visor	For 2-piece safety goggles, made of plastic	-	25 units
D	09498	Face visor	- Low reflection - With elastic strap	-	25 units
D	09483	Interlock full-face visor	- Low reflection - Distortion-free - With elastic strap - Latex-free - Fibreglass-free	-	24 units

### NEW ClearVizor® protective face shield set

Set comprising:

- 1 x frame
- 10 x foil shields
- Outstanding transparency
- No distortion of the field of vision
- Offers health protection when working
- Anti-glare
- Anti-static
- Equipped with anti-fogging properties
- Scratch-resistant
- Resistant to liquids and chemicals
- Safe (resistant to breaking and chipping)
- The visor is protected on both sides with adhesive foils (incl. removal aids)
- Light and ergonomic

	REF	Colour	SU
E	09714	White	1 set
F	09715	Blue	1 set
G	09716	Transparent	1 set
H	09717	Pink	1 set
I	09718	Yellow	1 set
J	09719	Green	1 set





A



B



C



D



E



F



G

**NEW ClearVizor® protective face shield set**

Set comprising:

- 1 x frame
- 2 x foil shields
- Outstanding transparency
- No distortion of the field of vision
- Offers health protection when working
- Anti-glare
- Anti-static
- Equipped with anti-fogging properties
- Scratch-resistant
- Resistant to liquids and chemicals
- Safe (resistant to breaking and chipping)
- The visor is protected on both sides with adhesive foils (incl. removal aids)
- Light and ergonomic

	REF	Colour	SU
A	09720	White	1 set
B	09721	Blue	1 set
C	09722	Transparent	1 set
D	09723	Pink	1 set
E	09724	Yellow	1 set
F	09725	Green	1 set

**NEW Accessories for ClearVizor®**

	REF	Description	SU
G	09726	Foil shield, without frame	10 units



A

B



C


D

**NEW InterClogs®**

- Made in one piece for residue-free reprocessing
- Anatomically shaped footbed for maximum comfort
- Anti-slip outsole for safety on smooth surfaces
- Material: TPE
- Heat resistance:
  - Water: 100°C (can be washed in a washer-disinfector)
  - Steam: 134°C
  - Dry heat: 120°C

The clogs can be fitted with different coloured straps.

	REF	Size	Colour	SU
A	26400	35/36	Navy blue	1 pair
A	26401	37/38	Navy blue	1 pair
A	26402	39/40	Navy blue	1 pair
A	26403	41/42	Navy blue	1 pair
A	26404	43/44	Navy blue	1 pair
A	26405	45/46	Navy blue	1 pair
B	26406	35/36	Green	1 pair
B	26407	37/38	Green	1 pair
B	26408	39/40	Green	1 pair
B	26409	41/42	Green	1 pair
B	26410	43/44	Green	1 pair
B	26411	45/46	Green	1 pair
C	26412	35/36	Red	1 pair
C	26413	37/38	Red	1 pair
C	26414	39/40	Red	1 pair
C	26415	41/42	Red	1 pair
C	26416	43/44	Red	1 pair
C	26417	45/46	Red	1 pair
D	26418	35/36	Yellow	1 pair
D	26419	37/38	Yellow	1 pair
D	26420	39/40	Yellow	1 pair
D	26421	41/42	Yellow	1 pair
D	26422	43/44	Yellow	1 pair
D	26423	45/46	Yellow	1 pair



A

B

D

C



E



F



G



H



I



J



K



L

**NEW InterClog® straps**

- Material: TPE
- Heat resistance:
  - Water: 100°C (can be washed in a washer-disinfector)
  - Steam: 134°C
  - Dry heat: 120°C

	REF	Size	Colour	SU
A	26424	35 - 40	Navy blue	1 pair
B	26425	35 - 40	Green	1 pair
C	26426	35 - 40	Red	1 pair
D	26427	35 - 40	Yellow	1 pair
A	26428	41 - 46	Navy blue	1 pair
B	26429	41 - 46	Green	1 pair
C	26430	41 - 46	Red	1 pair
D	26431	41 - 46	Yellow	1 pair

**Stainless steel dispenser for gloves**

	REF	Dimensions in mm	Designation	Description	SU
E	09200	262 x 152 x 90	Stainless steel dispenser for gloves	For tabletop or wall mounting, screws and plugs included	1 unit
F	09201	258 x 392 x 96	"TRIO" stainless steel dispenser for gloves	For mounting on the wall, screws and plugs included	1 unit
G	09202	250 x 135 x 95	"SOLO" stainless steel dispenser for gloves	For mounting on the wall, screws and plugs included	1 unit
H	20439	202 x 23 x 37	Standard rail adapter for glove dispenser	REF 09201 "Trio" and REF 09202 "Solo"	1 unit
I	09204	400 x 135 x 90	Stainless steel dispenser for gloves	For wall mounting, screws and plugs included, can be mounted vertically or horizontally, ideal for: REF 09588-09592	1 unit

**Shoe cover dispenser and accessories**

	REF	Dimensions in mm	Designation	Specification	SU
J	09965	455 x 235 x 165	Shoe cover dispenser	Aluminium housing	1 unit
K	09964	16 x 51	LDPE shoe covers for shoe cover dispenser – refill pack	<ul style="list-style-type: none"> <li>• 20 µ</li> <li>• 9 packs of 200</li> <li>• For single use</li> </ul>	1800 units

**Neck warmers with turtle neck**

- 100% cotton
- Anti-pilling
- Preshrunk (does not shrink when washed)
- Öko-Tex Standard 100

	REF	Size	SU
L	09439	S	1 unit
L	09403	M	1 unit
L	09404	L	1 unit




A



B



C



D



E



F



G

**Mask**

- 3-ply
- Blue and white
- With elastic ear loops
- For single use only

	REF	Visor	SU
A	09440	No	50/box
B	09441	Yes	20/box

**LDPE apron**


For single use only

	REF	Dimensions in mm	Description	SU
C	09495	70 x 136	- White - Material thickness: 50 µ	500 units

**Ear protectors**

	REF	Designation	Description	SU
D	09432	FLEXICAP ear protector	FLEXICAP headband ear protector with flexible ends, for single use only (replacement earplugs see REF 09434)	1 unit
D	09433	EARCAPS ear protector	EARCAPS headband ear protector with rigid band, for single use only (replacement earplugs see REF 09434)	1 unit
E	09434	EARCAPS replacement earplugs	Compatible with REF 09432 and 09433, for single use only	20 units
F	09435	EAR Classic 2 dispenser for earplugs	Simple to use: The earplugs are dispensed directly into your hand when the funnel-shaped container is twisted. With free access from the front or side, the container can be turned with one hand both clockwise and anti-clockwise. Incl. 1000 earplugs/dispenser	1 dispenser
G	09444	“EARSoft Yellow Neons” earplugs	Made of foam which slowly returns to its original shape For single use only, package: 250 pairs	250 pairs









**A** Kirschner wires made of stainless steel – trocar / round

REF	Length in mm	Dia. in mm	SU
07005	60	1.2	10 units
07010	60	1.8	10 units
07018	100	0.8	10 units
07020	100	1.0	10 units
07022	100	1.2	10 units
07403	150	0.7	10 units
07404	150	0.8	10 units
07406	150	1.0	10 units
07408	150	1.2	10 units
07409	150	1.4	10 units
07410	150	1.5	10 units
07411	150	1.6	10 units
07412	150	1.7	10 units
07413	150	1.8	10 units
07414	150	1.9	10 units
07415	150	2.0	10 units
07416	150	2.2	10 units
07417	150	2.4	10 units
07418	150	2.5	10 units
07035	160	0.8	10 units
07036	160	0.9	10 units
07037	160	1.0	10 units
07038	160	1.1	10 units
07039	160	1.2	10 units

REF	Length in mm	Dia. in mm	SU
07040	160	1.4	10 units
07042	160	1.6	10 units
07044	160	1.8	10 units
07046	160	2.0	10 units
07047	160	2.2	10 units
07048	160	2.4	10 units
07049	160	2.5	10 units
07051	160	3.0	10 units
07052	310	0.8	10 units
07053	310	0.9	10 units
07054	310	1.0	10 units
07055	310	1.1	10 units
07056	310	1.2	10 units
07057	310	1.4	10 units
07058	310	1.5	10 units
07059	310	1.6	10 units
07061	310	1.8	10 units
07063	310	2.0	10 units
07064	310	2.2	10 units
07065	310	2.4	10 units
07066	310	2.5	10 units
07067	310	2.8	10 units
07068	310	3.0	10 units



## A Kirschner wires made of stainless steel – trocar / trocar

REF	Length in mm	Dia. in mm	SU
07174	70	0.6	10 units
07175	70	0.8	10 units
07117	100	0.8	10 units
07118	100	0.9	10 units
07119	100	1.0	10 units
07121	100	1.2	10 units
07122	100	1.4	10 units
07172	150	0.7	10 units
07173	150	0.8	10 units
07134	160	0.8	10 units
07136	160	1.0	10 units
07138	160	1.2	10 units
07139	160	1.4	10 units
07141	160	1.6	10 units
07143	160	1.8	10 units
07147	160	2.4	10 units
07148	160	2.5	10 units
07150	160	3.0	10 units
07155	310	1.2	10 units
07156	310	1.4	10 units
07158	310	1.6	10 units
07159	310	1.7	10 units
07160	310	1.8	10 units
07162	310	2.0	10 units
07163	310	2.2	10 units



A



B



C



D



E



F



Kirschner wires made of stainless steel – trocar / flat

	REF	Length in mm	Dia. in mm	SU
A	07233	100	3.0	10 units
A	07234	160	0.8	10 units
A	07235	160	0.9	10 units
A	07236	160	1.0	10 units
A	07238	160	1.2	10 units
A	07239	160	1.4	10 units
A	07241	160	1.6	10 units
A	07242	160	1.7	10 units
A	07243	160	1.8	10 units
A	07245	160	2.0	10 units
A	07255	310	1.2	10 units
A	07257	310	1.5	10 units
A	07259	310	1.7	10 units
A	07260	310	1.8	10 units
A	07262	310	2.0	10 units
A	07263	310	2.2	10 units
A	07264	310	2.4	10 units
A	07265	310	2.5	10 units
A	07267	310	3.0	10 units

Miscellaneous Kirschner wires made of stainless steel

	REF	Length in mm	Dia. in mm	Ends	SU
B	07269	160	1.0	Trocar with 15 mm thread / round	10 units
C	07270	310	2.4	Drill tip / round	10 units
D	07268	430	2.4	Trocar with lengthwise knurl / special loop	10 units

Kirschner wires made of titanium – trocar / round

	REF	Length in mm	Dia. in mm	Colour	SU
E	07071	100	1.0	Silver	10 units
E	07072	100	1.2	Silver	10 units
E	07073	100	1.4	Silver	10 units
E	07074	100	1.6	Silver	10 units
E	07075	100	1.8	Silver	10 units
E	07076	100	2.0	Silver	10 units
E	07423	310	1.2	Silver	10 units
E	07422	310	1.8	Silver	10 units
F	07470	150	1.0	Gold	10 units
F	07471	150	1.2	Gold	10 units
F	07472	150	1.8	Gold	10 units





A

B

C

D

E



F

G

H

I

**Mesh sleeves for Kirschner wires**

Suitable for washer-disinfectors, heat-resistant

		REF	Length in mm	Dia./width in mm	Cross-section	Description	Material	SU
NEW	A	07473	170	20	Square	1 chamber	Stainless steel	1 unit
NEW	A	07474	200	20	Square	1 chamber	Stainless steel	1 unit
NEW	A	07475	330	20	Square	1 chamber	Stainless steel	1 unit
NEW	B	07476	170	20	Square	4 chambers Outlet opening: 0.8 - 1.0 - 1.2 - 1.4	Stainless steel	1 unit
NEW	B	07478	170	20	Square	4 chambers Outlet opening: 1.4 - 1.6 - 1.8 - 2.0	Stainless steel	1 unit
NEW	B	07477	200	20	Square	4 chambers Outlet opening: 0.8 - 1.0 - 1.2 - 1.4	Stainless steel	1 unit
NEW	B	07497	330	20	Square	4 chambers Outlet opening: 0.8 - 1.0 - 1.2 - 1.4	Stainless steel	1 unit
NEW	B	07499	330	20	Square	4 chambers Outlet opening: 1.2 - 1.4 - 1.6 - 1.8	Stainless steel	1 unit
NEW	B	07501	330	20	Square	4 chambers Outlet opening: 2.0 - 2.2 - 2.5 - without marking	Stainless steel	1 unit
	C	07430	170	22	Round	1 chamber, mesh size: 0.8 x 0.8 mm	Stainless steel	1 unit
	C	07431	200	22	Round	1 chamber, mesh size: 0.8 x 0.8 mm	Stainless steel	1 unit
	C	07461	330	22	Round	1 chamber, mesh size: 0.8 x 0.8 mm	Stainless steel	1 unit
	D	07460	170	20	Round	4 chambers, mesh size: 0.8 x 0.8 mm Outlet opening: 0.8 - 1.0 - 1.2 - 1.4	Stainless steel	1 unit
	E	07465	100	10	Round	1 chamber	Stainless steel	1 unit
	E	07400	160	10	Round	1 chamber	Stainless steel	1 unit
	E	07466	180	10	Round	1 chamber	Stainless steel	1 unit
	E	07467	200	10	Round	1 chamber	Stainless steel	1 unit
	E	07401	240	10	Round	1 chamber	Stainless steel	1 unit
	E	07468	250	10	Round	1 chamber	Stainless steel	1 unit
	E	07402	320	10	Round	1 chamber	Stainless steel	1 unit
	E	07421	440	10	Round	1 chamber	Stainless steel	1 unit
NEW	F	07493	130	20	Round	4 chambers, can be lasered as required	Titanium	1 unit
NEW	G	07494	200	20	Round	4 chambers, can be lasered as required	Titanium	1 unit
NEW	H	07495	350	20	Round	4 chambers, can be lasered as required	Titanium	1 unit
	I	07464	-	-	Round	Lid for REF 07430, 07431, 07461	Stainless steel	1 unit



A



B



C



D

**Kögel dispenser for Kirschner wires**

Made of stainless steel

	REF	Format in mm	For wire dia. in mm	For wire length in mm	Description	SU
A	08441	166-326 x 96 x 40	0.8 - 2.0	150 - 310	Flexible, incl. hex key, with plastic impact protection on both sides for 7 wires	1 unit
A	08439	166-326 x 96 x 40	1.4 - 3.0	150 - 310	Flexible, incl. hex key, with plastic impact protection on both sides for 7 wires	1 unit
B	08440	166 x 72 x 40	0.8 - 1.4	150	With plastic impact protection on both sides, for 5 wires	1 unit
C	08800	235 x 74 x 33	-	Up to 150	Set comprising: • 1 small parts sieve 1/6 • 1 lid for small parts sieve 1/6 • 1 wire dispenser	1 set

**Accessories for Kirschner wires**

	REF	Description	SU
D	07500	Gauge with wire diameters: 0.7 - 3.0 mm	1 unit



A



B

**Interlock washing sleeves made of stainless steel for Kirschner wires**

incl. 2 silicon caps, 1 ring holder and 1 tag – outside dia.: 19 mm – inside dia.: 17 mm

	REF	Length in mm	SU
A	07433	55	1 unit
A	07434	110	1 unit
A	07435	165	1 unit
A	07436	220	1 unit
A	07437	275	1 unit
A	07438	330	1 unit
A	07439	385	1 unit
A	07440	440	1 unit
A	07441	495	1 unit
A	07442	550	1 unit
A	07443	605	1 unit
A	07444	660	1 unit
A	07445	715	1 unit
A	07446	770	1 unit
A	07447	825	1 unit
A	07448	880	1 unit
A	07449	935	1 unit

**Interlock silicone cap for Interlock washing sleeves**

transparent – dia.: 19 mm – height: 7 mm

	REF	Description	SU
B	07451	Compatible with REF 07433 - 07449	10 units

<b>A</b>	
abox disposal containers	31
Accessories for ClearVizor	73
Accessories for Kirschner wires	91
Aluminium bowls	
Non-sterile	27
Sterile	27
<b>B</b>	
Bags	43
Beakers	29
Bed identification label	51
Bed labels	51
Boards made of POM synthetic material	51
Bowls & containers	21–31
Brush dispenser holders	55
<b>C</b>	
Cannulae	15–19
Carrying bag	43
C-bow cover	51
Clamp cover	45
Cleaning gloves	69
ClearVizor®	71–73
ClearVizor® protective face shield set	71–73
<b>D</b>	
Disposable bowls	
Plastic	23–27
Non-sterile	23, 27
Sterile	25, 27
Metal	27
Disposable cannulae	
Metal	19
Plastic	17
Disposable hand scrub brushes	55
Disposable pillow case	51
Disposable sharps containers	31
Disposal container	31
Dustproof bags	47
<b>E</b>	
EARCAPS	79
EAR Classic dispenser	79
Ear protectors	79
EARSoft	79
<b>F</b>	
Face visor	71
Fixing hand	49
Fixing hand accessories	49
FLEXICAP	79
Full-face visor	71
<b>G</b>	
Goggles frame	71
Gowns	59
Griffin beaker	29

<b>H</b>	
Hand scrub brushes	55
Hand scrub brushes and brush dispensers	53–55
Heat protective gloves	61
<b>I</b>	
Instrument trays	
Plastic	29
Stainless steel	29
InterClogs®	75
InterClog® straps	77
Interlock bags	43
Interlock disposable bowls	
Plastic	
Non-sterile	23
Sterile	25
Interlock disposable cannulae	
Metal	19
Plastic	17
Interlock full-face visor	71
Interlock goggles frame	71
Interlock gowns	59
Interlock MagnoMat	47
Interlock needle counters	47
Interlock PP sieve tray tag	39
Interlock repair tag	
Plastic	41
Tyvek	41
Interlock sharps pads	47
Interlock sieve tray tags	
Printed	37
Blank	35
Reel	39
Interlock spatulas	27
Interlock Surgical Scrub	55
Interlock suture bag	43
Interlock visor	71
Interlock washing sleeves	93
Interlock full-face visor	71
Interlock washing sleeves	93
<b>K</b>	
Kidney dishes	27
Kirschner wires	
Stainless steel	83–87
Titanium	87
Kirschner wires and accessories	81–93
<b>L</b>	
LDPE apron	79
LDPE shoe covers	77
Lid for mesh sleeves	89
Light handle covers	51
<b>M</b>	
MagnoMat	47
Mask	79
Medi scrub	55

## M

Mesh sleeves	
Round	89
Square	89
Metal bowls	
Non-sterile	27
Sterile	27
Microscope cover	51
Mini plastic cartridges	29

## N

Neck warmers	77
Needle counters	47

## O

OR and patient care	33–51
OR mat	47
OR microscope cover for Zeiss	51

## P

Pillow case	51
Plastic bowl	
Non-sterile	23
Sterile	25
PP sieve tray tags	39
PRObox disposal container	31
Protective face shields	71–73
Protective gloves	61–67
Mixture of neoprene and natural latex	63
Nitrile	65–67
Nomex	63
Protective sleeves	69

## R

Repair tags	
Plastic	41
Tyvek	41
RFID sieve tray tags	39

## S

Safety pins	47
Sharps pads	47
Shoe cover dispenser and accessories	77
Sieve tray tags	
Plastic	35–39
Blank	35
Reel	39
Special paper	35

## S

Silicone clamp cover	45
Silicone ring	49
Spatulas	27
Stainless steel dispenser for gloves	77
Stainless steel lid for mesh sleeves	89
Stainless steel rack for sieve tray tags	39
Surgical scrub	55
Suture bags	43

## T

Trocar / flat	87
Trocar / round	83, 87
Trocar / trocar	85

## V

Vessel loops	45
Visor	71

## W

Washing sleeves	93
Waste bags	43
Wire dispensers	91














**Interlock Medizintechnik GmbH**

Zum Windpark 1 • 23738 Lensahn • Germany

T: +49 4363 905900 • F: +49 4363 90590590

info@interlockmed.com • www.interlockmed.com

