

Our principles of Safety,  
Value and Performance make us  
the **No.1 choice worldwide**

**ACTIVE**  
BIRTHPOOLS  
[www.activebirthpools.com](http://www.activebirthpools.com)

# Active Birth Pools



  
**building  
better  
healthcare  
awards 2015**  
*Highly Commended*

The world's leading supplier  
of water birth pools - est. 1987

[www.activebirthpools.com](http://www.activebirthpools.com)



**Keith Brainin**  
Founder & Director

If you're looking for a water birth pool we're uniquely qualified to be of service

We helped pioneer the use of water for labour & birth by providing the pools that made this revolution in maternity care possible.

Since 1987 we've supplied 1,000's of water birth pools to hospitals world-wide.

Our approach is defined by a passion for excellence and a firm commitment to helping mothers & midwives

Mothers benefit from complete freedom of movement, unparalleled comfort, support and safety

Midwives report better results and higher levels of all-around satisfaction

We've achieved success by keeping things simple & focusing on our core principles of Safety, Value & Performance

### Our new pools are a quantum leap forward in every respect

1. Designed for optimum performance
2. Built to last for decades
3. Guaranteed for life
4. No.1 safety rating
5. Evidence based ergonomic design optimises experience
6. Maximum comfort & support for mothers and midwives.
7. Superior material enhances performance & durability



*Highly Commended*


*Highly commended for the quality of design and manufacture at the Building Better Healthcare Awards*

*The judges said that the Active Birth Pool succeed from all critical standpoints*





# The Active Birth Pool



Cad plan available

The pinnacle of water birth pool development.

An incomparable combination of superior materials & state-of-the-art manufacturing that provides mothers & midwives with unparalleled comfort, safety & support.

Specifications:	
<b>Dimensions:</b>	1840mm (l) x 1650mm (w) x 750mm (h) / 72" (l) x 65" (w) x 29.5" (h)
<b>Depth:</b>	640mm / 25"
<b>Weight (empty):</b>	90kg / 198lbs
<b>Capacity:</b>	550 litres (filled to 10cm below rim) / 145 gallons
<b>Material:</b>	Ficore® composite




[Click here](#) for further information

For more information please contact us or visit [www.activebirthpools.com](http://www.activebirthpools.com)



## The Princess Water Birth Pool



Cad plan available

Designed to fit in rooms where space is at a premium, the Princess Pool is surprisingly spacious and gives mothers room to explore and be supported in the most comfortable and beneficial positions.

### Specifications:

<b>Dimensions:</b>	1600mm (l) x 1100mm (w) x 770mm (h) / 63" (l) x 43" (w) x 30" (h)
<b>Depth:</b>	640mm / 25"
<b>Weight (empty):</b>	70kg / 154lbs
<b>Capacity:</b>	450 litres (filled to 10cm below rim) / 118 gallons
<b>Material:</b>	Ficore® composite



[Click here](#) for further information

For more information please contact us or visit [www.activebirthpools.com](http://www.activebirthpools.com)


5



# The Venus Water Birth Pool

The mid-size Venus Pool's design heritage stems from the proven principles we developed in the 1980's and takes them to new levels with seamless construction in Ficore® composite.

Specifications:	
Dimensions:	1920mm (l) x 1350mm (w) x 750mm (h) / 75.5" (l) x 65" (w) x 29.5" (h)
Depth:	640mm / 25"
Weight (empty):	80kg / 176lbs
Capacity:	500 litres (filled to 10cm below rim) / 132 gallons
Material:	Ficore® composite



Cad plan available



[Click here](#) for further information




# mothers & midwives love our water birth pools





*"Just writing to say thanks, I've worked with Active Birth Pools for many years. I always thought your pools were great but your new pools are amazing, make my job so much easier"*

G. RODINGHAM - NHS Midwife





For more information please contact us or visit [www.activebirthpools.com](http://www.activebirthpools.com)



# Light & Sound

## Light and sound play a significant role in the creation of a better birthing environment

Our easy-to-use lighting & sound systems make it simple for mothers & midwives to control and personalise the atmosphere & ambiance of the room.

In addition to making it easier for midwives to see into the pool the clear LED light illuminates the room with a soft glow.

The multi-colour LED lighting system allows users to cycle through a range of colours to find the hue most they find most appealing.

The bluetooth sound system enables users to wirelessly connect their phones or devices and listen to the music of their choice.


Two integral speakers turn the birth pool into a highly refined audio loudspeaker.



[Click here](#) for further information



# Features



For more information please contact us or visit [www.activebirthpools.com](http://www.activebirthpools.com)



## Step Unit

This sturdy step unit provides a safe, practical way for mothers to enter the pool

It gives the mother a height advantage.

She can step up – sit on the extra-wide rim and swivel into the water.

There is no climbing – she is grounded at all times & safeguarded from risk of slipping or falling.



[Click here](#) for further information

**Dimensions:** 455mm (l) x 300mm (w) x 150mm (h)  
/ 18" (l) x 12" (w) x 6" (h)

**Weight:** 3.85kg / 8.5lbs

Our principles of Safety,  
Value and Performance make us  
the **No.1 choice worldwide**

*"I can't tell you how much your pools have revolutionised birth."*

K. WILKIE - NHS Midwife

*"I've used the Active Birth pools for 15 years and I have always found them really beautiful, technically perfect & easy to use.*

*I love the simplicity of the lines that makes cleaning & disinfecting easy. I can't imagine a birth centre without an active birth pool!"*

DR. LUCIA ROCCA-IHENACHO - Lecturer in Midwifery & NIHR Research Fellow

*"We love your active birth pools, so much so we have one in every room."*

JO TALBOT - Chorley & Preston Birth Centre Manager

*"Most of the environments I see where women use water for labour & birth have Active Birth Pools.*

*Not only are they gloriously spacious & inviting, they are made to the highest specifications & considerations to maximise safety. In short, they are second to none!"*

SHEENA BYROM - OBE Midwife Speaker Author Honorary Fellow RCM Honorary Doctorate - Bournemouth University

TO BE A VIDEO



#### Active Birth Pools UK & International

+44 (0) 77 937 45907

keith@activebirthpools.com

#### Active Birth Pools America

651-895-2683

tim@activebirthpoolsusa.com

#### Active Birth Pools Australia/New Zealand

+61 (0) 2 9890 5687

cathy@pregnancy.com.au

#### Active Birth Pools Scandinavia

+45 4557 0881

bv@gm-medical.com

#### Active Birth Pools Benelux

+31 (0) 30 285 11 05

info@medicadynamics.nl

#### Active Birth Pools UAE

+971 (0) 4 259 1167

george@biomedic-intl.com

#### Active Birth Pools China

+400 1818 198

cnzcs66@163.com



# Active Birth Beanbag and Mat

The introduction of these few simple, easy to use, cost effective items transforms the birthing environment.

Specially designed to give mothers the freedom to move and be supported in the postures most beneficial for encouraging a physiological labour & natural, active birth.

The Active Birth Mat is perfect for use alongside a water birth pool.

It provides mothers with firm support for standing or squatting, cushions them when sitting or kneeling and is comfortable for lying down on.

The Active Birth Beanbag is taller & firmer than other beanbags to provide support for mothers in upright postures.

The medical grade waterproof fabric used to cover the mat & beanbag is soft and pleasing to the touch.



## Active Birth Beanbag

**Dimensions:** 800mm (dia.) x 700mm (h)  
/ 31.5" (dia) x 27.5" (h)

**Weight:** 4.48kg / 10lbs

## Active Birth Mat

**Dimensions:** 1500mm (l) x 1000mm (w) x 20mm (h)  
/ 59" (l) x 39.3" (w) x 0.75" (h)

**Weight:** 0.98kg / 2.2lbs



For more information please contact us or visit [www.activebirthpools.com](http://www.activebirthpools.com)