

PAJUNK®

InFiltr

*Cannula for the local
infiltration analgesia (LIA)*

NEW




Pain therapy

MADE IN GERMANY

InFiltr

Cannula for the local infiltration analgesia (LIA)

The local infiltration analgesia is performed by surgeons after a knee and hip joint placement (TEP) in a Single Shot procedure. As part of the fast-track concept the LIA is used for early rehabilitation. With InFiltr PAJUNK® offers an infiltration cannula with a large infiltration field. Due to the closed tip and the six lateral openings the anaesthetic can be widely dispersed in the surgical wound.



Advantages

- ➔ *Easy puncture*
- ➔ *Greater infiltration field*
- ➔ *Closed tip and six lateral openings*
- ➔ *Early mobilisation*
- ➔ *Shorter hospitalisation*

Product	Size	Item no.	PU
InFiltr Infiltration cannula with closed tip and six lateral openings on the first 20 mm	20G x 100 mm	1158-3G100	25

PAJUNK® GmbH
Medizintechnologie
Karl-Hall-Strasse 1
D-78187 Geisingen / Germany
Phone +49 (0) 77 04/92 91-0
Telefax +49 (0) 77 04/92 91-6 00
www.pajunk.com

PAJUNK® Medical Produkte GmbH
D.A.CH • BeNeLux
Karl-Hall-Strasse 1
D-78187 Geisingen / Germany
Phone +49 (0) 77 04/80 08-0
Telefax +49 (0) 77 04/80 08-150
www.pajunk.com