

perfect quality · sincere service



Tappa 坦帕

杭州坦帕医疗科技有限公司
Hangzhou Tappa Medical Technology Co., Ltd.



Hangzhou Tappa Medical Technology Co., Ltd.

Address: NO. 225#, Chutian Road, BinJiang District, Hangzhou, China, 310051.

Office: (+86) 571 8808 5616

Fax: (+86) 571 8808 5016

Website: www.tappamed.com



ABOUT US

Hangzhou Tappa Medical Technology Co., Ltd. was established in 2008. The headquarter and production & R&D center are located in Binjiang district, Hangzhou city. The company is mainly specialized in research, production, and sales of anesthesia consumables in anesthesia, breathing, ICU, operating theatre, etc. The products are widely used in anesthesia, therapy and emergency.

Tappa Medical is one of the main anesthesia consumable suppliers in the world. Tappa had set up business relationship with hundreds of medical device distributors all over the world and a lot of domestic exporters. The products are already sold to Europe, Middle East, Asia, Africa which have big market share and more strict requirements. In domestic market, our products are sold to all over the country and used in over 400 First Class Top Level hospitals.

Now Tappa had passed ISO 13485 & CE, and FDA is on the way.

Tappa is in a new period of development. Adhering to the "perfect quality, sincere service" purposes, Tappa will consistently support the development of medical industry, and strive for technological innovation, academic support, nurturing talent and other aspects.

CATALOG

01 Laryngeal Mask Airway

(Single Use & Reusable)

- Two Way Laryngeal Mask Airway
- Classic One Way Silicone Laryngeal Mask Airway
- One Way Steel Reinforced Laryngeal Mask Airway

02 Endotracheal Tube

- Silicone Reinforced Endotracheal Tube
- Bite Block Silicone Reinforced Endotracheal Tube
- PVC Wire Reinforced Endotracheal Tube
- PVC Endotracheal Tube
- PVC Oral Endotracheal Tube
- PVC Nasal Endotracheal Tube
- Endotracheal Tube with Evacuation Lumen

03 Double Lumen Endobronchial Tube

- Silicone Double Lumen Endobronchial Tube
- PVC Double Lumen Endobronchial Tube

04 Endobronchial Blocker Tube

- Auto-Inflated Endobronchial Blocker Tube
- Endobronchial Blocker Tube



05 Breathing Circuits

- Breathing Circuit
- Anesthesia Circuit

06 Catheter Mount

07 Breathing System Filter

- Bacterial Virus Filter
- HME Filter

08 Nasopharyngeal Airway

09 Oropharyngeal Airway

10 Disposable Air Cushion Face Mask


11 Tracheal Tubes Holder

12 Intubating Stylet

13 Tracheal Tube Kit

01 Laryngeal Mask Airway

(Single Use & Reusable)



Two Way Laryngeal Mask Airway

Two Way Laryngeal Mask Airway

- The built-in drain tube that allows expelled gastric content to bypass the pharynx. This specific feature is designed to decrease the risk of aspiration
- Made of medical silicone material, non-toxic & non-stimulating
- Steel reinforced airway tube can prevent biting & kinking
- Double layer cuff seals the throat perfectly
- Good biocompatibility
- Minimal haemodynamic response
- Reusable type can be used up to 40 times
- Package: 10pcs/box, 60pcs/carton

Size	Suitable Weight	Reference NO#
1.0	<5kg	LMA1004P2
1.5	5-10kg	LMA1507P2
2.0	10-20kg	LMA2010P2
2.5	20-30kg	LMA2514P2
3.0	30-50kg	LMA3020P2
4.0	50-70kg	LMA4030P2
5.0	>70kg	LMA5040P2

Classic One Way Silicone Laryngeal Mask Airway



Classic One Way Silicone Laryngeal Mask Airway

- Made of medical silicone material, non-toxic & non-stimulating
- Good biocompatibility
- Soft silicone cuff fits tightly with the throat
- Facilitates smooth emergence from anaesthesia
- Minimal haemodynamic response
- Reusable type can be used up to 40 times
- Package: 10pcs/box, 60pcs/carton

Size	Suitable Weight	Reference NO#
1.0	<5kg	LMA1004C2
1.5	5-10kg	LMA1507C2
2.0	10-20kg	LMA2010C2
2.5	20-30kg	LMA2514C2
3.0	30-50kg	LMA3020C2
4.0	50-70kg	LMA4030C2
5.0	>70kg	LMA5040C2

One Way Steel Reinforced Laryngeal Mask Airway

One Way Steel Reinforced Laryngeal Mask Airway



- Made of medical silicone material, non-toxic & non-stimulating
- Steel reinforced airway tube can prevent biting & kinking
- Good biocompatibility
- Soft silicone cuff fits tightly with the throat
- Facilitates smooth emergence from anaesthesia
- Minimal haemodynamic response
- Reusable type can be used up to 40 times
- Package: 10pcs/box, 60pcs/carton

Size	Suitable Weight	Reference NO#
1.0	<5kg	LMA1004F2
1.5	5-10kg	LMA1507F2
2.0	10-20kg	LMA2010F2
2.5	20-30kg	LMA2514F2
3.0	30-50kg	LMA3020F2
4.0	50-70kg	LMA4030F2
5.0	>70kg	LMA5040F2

02 Endotracheal Tube



Cuffed

Silicone Reinforced Endotracheal Tube (Cuffed & Uncuffed)

- Made of medical silicone material, non-toxic & non-stimulating, good biocompatibility
- Steel reinforced tube can prevent biting & kinking
- High volume low pressure cuff reduces the pressure to tracheal wall
- Accompanied with intubating stylet
- Package: 10pcs/box, 100pcs/carton



Uncuffed

Cuffed

Size	Reference NO#
3.5	ETT3524C
4.0	ETT4024C
4.5	ETT4524C
5.0	ETT5024C
5.5	ETT5524C
6.0	ETT6024C
6.5	ETT6524C
7.0	ETT7024C
7.5	ETT7524C
8.0	ETT8024C
8.5	ETT8524C
9.0	ETT9024C
9.5	ETT9524C

Uncuffed

Size	Reference NO#
2.5	ETT2524
3.0	ETT3024
3.5	ETT3524
4.0	ETT4024
4.5	ETT4524
5.0	ETT5024
5.5	ETT5524
6.0	ETT6024
6.5	ETT6524
7.0	ETT7024

Cuffed

Size	Reference NO#
3.5	ETT3525C
4.0	ETT4025C
4.5	ETT4525C
5.0	ETT5025C
5.5	ETT5525C
6.0	ETT6025C
6.5	ETT6525C
7.0	ETT7025C
7.5	ETT7525C
8.0	ETT8025C
8.5	ETT8525C
9.0	ETT9025C
9.5	ETT9525C

Uncuffed


Size	Reference NO#
3.5	ETT3525
4.0	ETT4025
4.5	ETT4525
5.0	ETT5025
5.5	ETT5525
6.0	ETT6025
6.5	ETT6525
7.0	ETT7025



Cuffed

Bite Block Silicone Reinforced Endotracheal Tube (Cuffed & Uncuffed)

- Made of medical silicone material, non-toxic & non-stimulating, good biocompatibility
- Bite block design can protect the teeth
- Steel wire reinforced tube can prevent biting & kinking
- High volume low pressure cuff reduces the pressure to trachea
- Accompanied with intubating stylet
- Package: 10pcs/box, 100pcs/carton



Uncuffed



Cuffed

E^{PVC} Wire Reinforced Endotracheal Tube (Cuffed & Uncuffed)

- Medical PVC, DEHP free
- Steel wire reinforced tube can prevent biting & kinking
- High volume low pressure cuff reduces the pressure to trachea
- Radio opaque line for x-ray visualization
- Murphy eye
- Package: 10pcs/box, 100pcs/carton




Uncuffed

Cuffed	
Size	Reference NO#
3.0	ETT3014C
3.5	ETT3514C
4.0	ETT4014C
4.5	ETT4514C
5.0	ETT5014C
5.5	ETT5514C
6.0	ETT6014C
6.5	ETT6514C
7.0	ETT7014C
7.5	ETT7514C
8.0	ETT8014C
8.5	ETT8514C
9.0	ETT9014C
9.5	ETT9514C
10.0	ETT10014C

Uncuffed	
Size	Reference NO#
3.0	ETT3014
3.5	ETT3514
4.0	ETT4014
4.5	ETT4514
5.0	ETT5014
5.5	ETT5514
6.0	ETT6014
6.5	ETT6514
7.0	ETT7014
7.5	ETT7514
8.0	ETT8014
8.5	ETT8514
9.0	ETT9014

Cuffed	
Size	Reference NO#
3.0	ETT3011C
3.5	ETT3511C
4.0	ETT4011C
4.5	ETT4511C
5.0	ETT5011C
5.5	ETT5511C
6.0	ETT6011C
6.5	ETT6511C
7.0	ETT7011C
7.5	ETT7511C
8.0	ETT8011C
8.5	ETT8511C
9.0	ETT9011C
9.5	ETT9511C
10.0	ETT10011C


Uncuffed	
Size	Reference NO#
2.0	ETT2011
2.5	ETT2511
3.0	ETT3011
3.5	ETT3511
4.0	ETT4011
4.5	ETT4511
5.0	ETT5011
5.5	ETT5511
6.0	ETT6011
6.5	ETT6511
7.0	ETT7011
7.5	ETT7511
8.0	ETT8011
8.5	ETT8511
9.0	ETT9011
9.5	ETT9511
10.0	ETT10011



Cuffed

E^{PVC} Endotracheal Tube (Cuffed & Uncuffed)

- Medical PVC, DEHP free
- Steel wire reinforced tube can prevent biting & kinking
- High volume low pressure cuff reduces the pressure to trachea
- Radio opaque line for x-ray visualization
- Murphy eye
- Package: 10pcs/box, 100pcs/carton



Uncuffed



Cuffed

PVC Oral Endotracheal Tube (Cuffed & Uncuffed)

- Medical PVC, DEHP free, transparent
- Oral preformed to reduce tube kinking during oral intubation
- Used for nasal surgery, ophthalmic surgery, facial surgery
- High volume low pressure cuff reduces the pressure to trachea
- Radio opaque line for x-ray visualization
- Murphy eye
- Package: 10pcs/box, 100pcs/carton



Uncuffed

Cuffed	
Size	Reference NO#
4.0	ETT4012C
4.5	ETT4512C
5.0	ETT5012C
5.5	ETT5512C
6.0	ETT6012C
6.5	ETT6512C
7.0	ETT7012C
7.5	ETT7512C
8.0	ETT8012C
8.5	ETT8512C
9.0	ETT9012C

Cuffed	
Size	Reference NO#
4.0	ETT4013C
4.5	ETT4513C
5.0	ETT5013C
5.5	ETT5513C
6.0	ETT6013C
6.5	ETT6513C
7.0	ETT7013C
7.5	ETT7513C
8.0	ETT8013C
8.5	ETT8513C
9.0	ETT9013C

Uncuffed	
Size	Reference NO#
3.0	ETT3013
3.5	ETT3513
4.0	ETT4013
4.5	ETT4513
5.0	ETT5013
5.5	ETT5513
6.0	ETT6013
6.5	ETT6513
7.0	ETT7013

Uncuffed	
Size	Reference NO#
3.0	ETT3012
3.5	ETT3512
4.0	ETT4012
4.5	ETT4512
5.0	ETT5012
5.5	ETT5512
6.0	ETT6012
6.5	ETT6512
7.0	ETT7012

Size	Reference NO#
6.0	ETT6014E
6.5	ETT6514E
7.0	ETT7014E
7.5	ETT7514E
8.0	ETT8014E
8.5	ETT8514E
9.0	ETT9014E



Cuffed

PVC Nasal Endotracheal Tube (Cuffed & Uncuffed)

- Medical PVC, DEHP free, transparent
- Nasal preformed to reduce tube kinking during nasal intubation
- Used for oral surgery, maxillofacial surgery
- High volume low pressure cuff reduces the pressure to trachea
- Radio opaque line for x-ray visualization
- Murphy eye
- Package: 10pcs/box, 100pcs/carton




Uncuffed

W Endotracheal Tube With Evacuation Lumen

- No trauma to vocal cords when suction
- Help to reduce the rate of VAP (Ventilator Associated Pneumonia)
- Medical grade PVC, Non-toxic, Latex free
- X- ray through the tube
- Smooth Murphy Eyelet
- High volume low pressure bolloon
- Package: 10pcs/box, 100pcs/carton



03 Double Lumen Endobronchial Tube



ESilicone Double Lumen Endobronchial Tube

- Made of medical silicone material, non-toxic & non-stimulating, good biocompatibility
- Used for one-lung ventilation
- Steel reinforced top can prevent kinking and movement, convenient for radiography location
- High volume low pressure cuff reduces the pressure to trachea
- Accompanied with intubating stylet, suction tube & connectors
- Package: 1pcs/box, 10pcs/carton

Size	Reference NO#
32Fr	DLT0132S
35Fr	DLT0135S
37Fr	DLT0137S
39Fr	DLT0139S



EPVC Double Lumen Endobronchial Tube (Left side & Right side & Left side with Carina)


- Made of Medical PVC, DEHP free, transparent
- Used for one-lung ventilation
- High volume low pressure cuff reduces the pressure to trachea
- Accompanied with intubating stylet, suction tube & connectors
- Package: 1pcs/box, 10pcs/carton

Size	Reference NO#(Left Side)	Reference NO#(Right Side)	Reference NO#(Left Side with Carina)
26Fr	DLT0126	DLT0226	DLT0326
28Fr	DLT0128	DLT0228	DLT0328
32Fr	DLT0132	DLT0232	DLT0332
35Fr	DLT0135	DLT0235	DLT0335
37Fr	DLT0137	DLT0237	DLT0337
39Fr	DLT0139	DLT0239	DLT0339
41Fr	DLT0141	DLT0241	DLT0341

04 Endobronchial Blocker Tube

Auto-inflation Endobronchial Blocker Tube


- Unique automatic inflation system enables the operator to inflate the cuff with one hand and operate bronchi fiberscope simultaneously
- Endobronchial Blocker Set allows one-lung ventilation by using a conventional endotracheal tube and a fiberoptic bronchoscope
- No need to replace trachea tubes
- The stability of cuff balloon's internal pressure is better than any other products, and the low pressure cuffs minimize trauma
- Package: 1 set/box, 10 sets/carton



Size	Reference NO#
7Fr	EBT0207
9Fr	EBT0209

Endobronchial Blocker Tube (Without Auto-inflation)

- Endobronchial Blocker Tube allows one-lung ventilation by using a conventional endotracheal tube and a fiberoptic bronchoscope
- No need to replace trachea tubes
- The stability of cuff balloon's internal pressure is better than any other products, and the low pressure cuffs minimize trauma
- Package: 1 set/box, 10 sets/carton



Size	Reference NO#
5Fr	EBT0105
7Fr	EBT0107
9Fr	EBT0109

05 Breathing Circuits

Circuits for breathing machine

(with water traps, pediatric & adult)

Circuits for anesthesia machine

(1.2m, 1.5m, 1.8m, expandable & unexpandable, pediatric & adult)



- All of our circuits are made with only premium-quality materials and have proven reliable during use

Type	Reference NO#	Original Length	Expandable Length
Anesthesia Circuits(M)	BTM0112E	0.43m	1.2m
Anesthesia Circuits(M)	BTM0115E	0.54m	1.5m
Anesthesia Circuits(M)	BTM0115C	0.54m	1.5m
Anesthesia Circuits(M)	BTM0118C	0.65m	1.8m
Anesthesia Circuits(M)	BTM0212E	1.2m	/
Anesthesia Circuits(M)	BTM0215E	1.5m	/
Anesthesia Circuits(M)	BTM0215C	1.5m	/
Anesthesia Circuits(M)	BTM0218C	1.8m	/
Breathing Circuits(H)	BTH0212E	1.2m	/
Breathing Circuits(H)	BTH0215E	1.5m	/
Breathing Circuits(H)	BTH0215C	1.5m	/
Breathing Circuits(H)	BTH0218C	1.8m	/

06 Catheter Mount

- Disposable Catheter Mount Double Swivel Elbow
- Including standard and collapsible tube
- Providing convenience for surgery or other clinical process
- A variety of connectors for airway management options



07 Breathing System Filter



Bacterial Virus Filter (BVF)

- Bacterial virus filter
- Breathing filters help eliminate the risk of cross-infection by preventing the transfer of bacteria and viruses to and from the patient
- Protect the patients and ventilators from cross contamination
- Low respiratory resistance & low dead space
- Package: 1pcs/pouch, 100pcs/Carton

Reference NO#	Capacity	Connector	Moisturization
HME101	92ml	22M/15F-22F/15M	No
HME102	92ml	22M-22F	No
HME103	62ml	22M/15F-22F/15M	No
HME104	62ml	22M/15F-22F	No
HME105	43ml	22M/15F-22F/15M	No
HME106	43ml	22M/15F-22F/15M	No

HMEF (Heat And Moisture Exchanger Filter)

- Good humidification function
- Low respiratory resistance & low dead space
- Types for pediatric, child, and adult
- Breathing filters help eliminate the risk of cross-infection by preventing the transfer of bacteria and viruses to and from the patient
- Protect the patients and ventilators from cross contamination
- Package: 1pcs/pouch, 100pcs/Carton



Reference NO#	Capacity	Connector	Moisturization
HME201	95ml	22M/15F-22F/15M	Yes
HME202	95ml	22M/15F-22F/15M	Yes
HME203	45ml	22M/15F-22F/15M	Yes
HME204	45ml	22M/15F-22F/15M	Yes
HME205	45ml	22M/15F-22F/15M	Yes
HME206	66ml	22M/15F-22F/15M	Yes
HME207	92ml	22M/15F-22F/15M	Yes
HME208	92ml	22M/15F-22F/15M	Yes
HME209	26ml	22M/15F-22F/15M	Yes
HME210	26ml	22M/15F-22F/15M	Yes
HME211	99ml	22M/15F-22F/15M	Yes
HME212	99ml	22M/15F-22F	Yes
HME213	5ml	22M/15F-22F/15M	Yes
HME214	5ml	22M/15F-22F	Yes
HME215	5ml	15F-15M	Yes
HME216	95ml	22M/15F-22F/15M	Yes
HME217	66ml	22M/15F-22F/15M	Yes

08 Nasopharyngeal Airway

- Made of medical PVC material, soft & comfortable
- Only for nasal airway



Size

- 2.5mm
- 3.0mm
- 3.5mm
- 4.0mm
- 4.5mm
- 5.0mm
- 5.5mm
- 6.0mm
- 6.5mm
- 7.0mm
- 7.5mm
- 8.0mm
- 8.5mm
- 9.0mm

09 Oropharyngeal Airway

- Medical polyethylene material
- color-coded by size for easy identification



Size

- 3.0cm
- 3.5cm
- 4.0cm
- 4.5cm
- 5.0cm
- 5.5cm
- 6.0cm
- 6.5cm
- 7.0cm
- 8.0cm
- 9.0cm
- 10.0cm
- 11.0cm
- 12.0cm

10 / Disposable Air Cushion Face Mask



- Surgical Disposable Anesthesia Breathing Mask
- The transparent mask shell help the medical care personel to see complicating disease like vomit as soon as possible
- Soft air cushion can keep mask sealed; Use only once, avoid cross-infection


Size
0 #
1 #
2 #
3 #
4 #
5 #
6 #

11 / Tracheal Tubes Holder



- Made of nontoxic medical-grade PP material
- Fits various sizes of endotracheal tubes
- Accomodate all types of laryngeal masks
- Full padded foam to back for patient's comfort
- Allows for in-use suctioning of the oropharynx

12 / Intubating Stylet



- Made of Malleable aluminium which allows the stylet to be easily adapted to any desired shape
- Soft Rounded tip for minimizing patient trauma during intubation
- Flexible use for assitance of intubation

Size
6Fr
10Fr
14Fr

13 / Tracheal Tube Kit



- Including endotracheal tube, suction tube, suction connecting tube, injection syringe, endotracheal tube holder, oropharyngeal airway, laryngoscope, gloves, intubating stylet and blanket
- Full products for use
- Different endotracheal tubes for choosing
- 1set/box, 20boxes/carton

REF NO	Type
EB-P	PVC Endotracheal Tube Kit
EB-S1	Reinforced PVC Endotracheal Tube Kit
EB-S2	Reinforced Silicone Endotracheal Tube Kit