



**PROMETHEUS
MEDICAL LTD**

**SPLINTAGE &
IMMOBILISATION**



www.prometheusmed.com



Prometheus Traction Splint

The Prometheus Traction Splint is designed to be used on those patients with a suspected limb fracture. Proper use can decrease the pain and damage caused by the fracture.

Key Features:

- One person application in less than 2 minutes
- Traction pole adjusts to different sizes instantly for adult or paediatric use
- No need to roll patient
- No uncomfortable Ischial bar
- Can be used with the Prometheus Pelvic Splint
- Lightweight and compact

Product Code:

PDF101



Prometheus Omni Splint®

The Prometheus Omni Splint® is a malleable closed cell foam splint with an aluminium backbone.

Key Features:

- Waterproof
- Latex free
- Lightweight and compact
- Reusable
- Secured in place with tape or bandage of choice
- Closed cell foam allows for simple cleaning and disinfection

Specification:

- Pack Size: Single
- Dimensions: 91cm x 11cm (open)
11.5cm x 8 cm (packed)

Product Code:
PDF105





Prometheus Pelvic Splint

Reduction and stabilisation of pelvic injuries should occur as soon as possible, on scene. Studies of patients sustaining pelvic injuries have shown that rapid immobilisation reduces haemorrhage by stabilising the fractured bone ends, reducing bleeding and stabilising clot formation.

The Prometheus Pelvic Splint offers a 'one size fits all' solution, which is light, easy to apply and highly effective.

Key features:

- Neoprene construction - very light
- One size fits all - from child to bariatric patient
- May be cut to size to allow access to abdomen and groin without compromise
- Soft material with no stitching reduces risk of pressure ulceration or friction
- Perineum and anus are left uncovered

Specification:

- Material - 6mm stretch Neoprene fabric bonded with loop nylon on both faces
- Neoprene Sling - 126cm x 19cm
- Hook and Loop Anchor - 10cm x 18cm
- Hook and Loop Belt - 43cm x 5cm
- Hook and Loop Tab - 10cm x 13cm
- Tri-glide Plastic Loop - 5cm

Product Code: PRO400

NSN: 6515-99-174-5924



Prometheus Medical Ltd The Old Rectory,
Hope Under Dinmore, Herefordshire HR6 0PW

Tel: +44 (0) 1568 613942

Fax: +44 (0) 1568 620032

www.prometheusmed.com