



● ● ●
TruPICC

Product Code: TPIC100

TruPICC weight: 6kg

TruPICC dimensions: 90 x 35 x 20cm

Full shipment weight: 8kg

Full shipment dimensions: 95 x 37 x 21cm

*Part of the TruUltra product range, **TruPICC** is an anatomical adult male arm designed for teaching the skills associated with ultrasound guided PICC line and Intravenous (IV) techniques.*

TruPICC develops the user's skills associated with needle placement, guidewires and catheters. The upper arm anatomy features the brachial, basilic, cephalic and median cubital veins, the superior vena cava and the brachial artery for vessel differentiation.

Ideal for training in...

- **Ultrasound guided PICC Line Insertion**

- Full catheter PICC Line placement and seldinger techniques
- Realistic visualisation of the vascular anatomy of an arm including the brachial artery, and the brachial, cephalic, basilic and median cubital veins
- Upper chest features the superior vena cava to verify correct catheter placement using ultrasound

- **Intravenous (IV) Cannulation Techniques**

- Realistic visualisation of the median cubital, brachial and basilic veins
- Differentiation between the basilic vein and the basilic artery
- Positive user feedback upon successful entry into the vessels by withdrawal of fluid
- Lifelike vascular 'tenting' upon successful entry into the vessel

Features and benefits...

- **Realistic:** Anatomically accurate vascular features including the brachial, cephalic, basilic and median cubital veins, and the brachial artery
- **Lifelike training:** Lifelike vascular 'tenting' upon successful entry into the vessels
- **Successful procedure verification:** Upper chest features superior vena cava to verify correct PICC catheter placement using ultrasound. Withdrawal of blood verifies correct IV cannulation
- **Quick Set-Up:** Initial set-up time is less than 5 mins
- **Durable:** Unique, self-regenerating TruUltra materials diminish the appearance of needle tracks to facilitate a high volume of repeated practice. The insert will facilitate approx. 800 needle penetrations
- **Interchangeable:** Replacement insert is cost-effective and can be replaced quickly when required
- **Transportable:** Delivered in a durable carry case for easy transportation and safe storage
- **Universal:** Compatible with all Ultrasound machine brands

Package Contents...

- 1 TruPICC model
- 1 TruCorp carry case
- 1 USB pen drive containing product user manual
- 1 bottle of artificial blood concentrate (250ml)
- 1 luerlock syringe for fluid system preparation/withdrawal

