# **ORDERING INFORMATION**

NO.	SIZE	QTY.
400150	6F	Box (10 each)
400120	7F	Box (10 each)
400160	8F	Box (10 each)
400180	10F	Box (10 each)
400200	12F	Box (10 each)

#### Additional SecurAcath<sup>®</sup> product information

- Not made with natural latex rubber
- MRI Conditional

More information available at www.securacath.com

#### **TO ORDER:**

Call customer service at 866-980-1811 or email order@securacath.com






X

security

HOID





181 Cheshire Lane, Suite 100 Plymouth, MN 55441 USA +1.763.225.6699 www.securacath.com


Please refer to instructions for use for indications, contraindications, hazards, warnings, cautions and directions for use.

CE 0413 www.securacath.com/patents ©2020 by Interrad Medical, Inc. Interrad Medical and SecurAcath are trademarks of Interrad Medical, Inc. All rights reserved.

1329-036 Rev C

# A New Securement Device for Drains

# SecurAcath<sup>®</sup> provides improved drainage catheter securement





- Small, blunt, nitinol securement feet are placed just beneath skin right at the catheter insertion site
- Cover is snapped onto base to affix to catheter shaft
- No sutures or additional skin punctures are needed
- No adhesives needed for securement
- Remains in place for life of catheter



FOLD



INSERT



**SNAP** 

## **Benefits on Drainage Catheters**

#### **Decreased Catheter Dislodgments** and Replacement Costs

- Lower cost of patient care
- Minimize therapy interruption

#### **Increased Stability**

- SecurAcath secures right at the catheter/skin junction
- Minimizes catheter movement, pistoning and migration

#### Eliminates Suture Needle Sticks

• No sharps related injuries

#### Improved Efficiency

SecurAcath is much faster to use than sutures

#### **360** Degree Site Cleaning While Secured

- Excellent cleaning access around the entire insertion site
- Catheter remains stable and secure during cleaning

#### **Decreased Skin Irritation**

- Sutures can become infected, cause skin irritation and breakdown
- Minimize use of skin-damaging adhesives

### SecurAcath<sup>®</sup> – For the Life of the Line

Find out more at www.securacath.com



