



FOR EVERY INTUBATION.  
EVEN THE STRANGE ONES.

The new HEINE visionPRO video laryngoscope.



# PRODUCT OVERVIEW

The new HEINE visionPRO video laryngoscope



# The new HEINE visionPRO video laryngoscope.

With very, very, very good image quality.

- ▶ Outstanding image quality
- ▶ Robust and smart product design
- ▶ Ecological & Economical





# SINGLE-USE BLADE MAC



## First-ever blade from upcycled material

- ▶ Upcycled material from discarded refrigerators
- ▶ New HEINE FixLock System
- ▶ Safe & hygienic
- ▶ Established Macintosh form in sizes 3 & 4



# CHARGING STATION



## Charge and store

- ▶ Ecological & Economical
  - Rechargeable battery
  - Recurring costs spent on buying disposable batteries are eliminated
- ▶ Smart charging management



# FEATURES & BENEFITS

The new HEINE visionPRO video laryngoscope

# OUTSTANDING IMAGE QUALITY



## Bright and sharp image

- ▶ Unique worldwide: HEINE transfective allBRIGHT display adapts to even the brightest ambient light
- ▶ Excellent 3.2" display in practical portrait format
- ▶ Low-reflection intubation thanks to an adjustable display angle
- ▶ Realistic color representation thanks to smart image processing
- ▶ HEINE LED<sup>HQ</sup>

## Quick and comprehensive orientation

- ▶ Situational field of view of the relevant glottis and the surrounding anatomy



# ROBUST AND SMART PRODUCT DESIGN

## Extremely durable and long-lasting

- ▶ Grip head and display frame made of anodised aluminium
- ▶ Stainless steel handle and camera arm
- ▶ Thick screen protector

## Reliable, accurate and intuitive operation – even in time-critical situations

## Safe: Display and handle can be separated for numerous validated reprocessing methods

- ▶ Handle: from wipe disinfection to immersion and low temperature sterilisation
- ▶ Display: low level / wipe disinfection reprocessing





# ECOLOGICAL

## Blade made from recycled material

- ▶ Use of plastics from recycled refrigerators
- ▶ Material savings through minimalist design
- ▶ Tons of savings in CO<sub>2</sub>

## Rechargeable battery vs. disposable batteries

- ▶ Reduced carbon footprint by avoiding the consumption of resources when manufacturing new batteries
- ▶ Reducing the consumption of rare and heavy metals
- ▶ Significant reduction in hazardous waste volume



# ECONOMICAL




## Cost reduction through economic procurement management

- ▶ Depending on your needs, displays and handles can be purchased separately
- ▶ Elimination of the recurring costs for disposable batteries





**MORE INFORMATION:  
[WWW.HEINE.COM](http://WWW.HEINE.COM)**



**HEINE** QUALITY  
MADE IN GERMANY

HEINE Optotechnik GmbH & Co. KG  
Dornierstr. 6 | 82205 Gilching | Germany  
Tel. +49 (0)8152 38-0 | [info@heine.com](mailto:info@heine.com) | [heine.com](http://heine.com)